

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | GEMFI<br>600MG,STKF2A00 |
| User ID :        | sj97904       | Batch NO. :          | KFM2500260              |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 13 Min 21 Sec   |
| Total Counter :  | 400070        | Total Running Time : | 3296 Hour 45 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 198               | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 600000            | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50                | Min     | Supply Time ;      | 0.6          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 10.21             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 9.82              | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 8.01              | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 7.86              | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 7.70              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 5.89              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 5.50              | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |                |       |            |      |       |            |    |
|----------------------|----------------|-------|------------|------|-------|------------|----|
| Mean Pre Pressure :  | 1.56           | kN    | Min :      | 1.40 | kN    | Punch NO.; | 15 |
| SD :                 | 6.20           | %     | Max :      | 1.90 | kN    | Punch NO.; | 22 |
| Mean Main Pressure : | 8.09           | kN    | Min :      | 7.71 | kN    | Punch NO.; | 52 |
| SD :                 | 2.30           | %     | Max :      | 8.60 | kN    | Punch NO.; | 17 |
| Filling Depth :      | Increase ; 317 | Times | Decrease ; | 312  | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 30       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.90    | mm  | Pre Thickness ;  | 3.91     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.55    | mm  | Main Thickness ; | 3.55     | mm           |
| Filling Depth :               | 8.31    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.775   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 18      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 19      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 44 |
| Lower Punch :    | Average ; | 0.28 | kN |             |    |
|                  | Max ;     | 1.04 | kN | Max Punch ; | 49 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 400070            | Tablets |       |            |         |
| A.W.C Counter : | 391290            | Tablets |       |            |         |
| HEP :           | 611 (0.2%)        | Tablets | LEP : | 510 (0.1%) | Tablets |
| GOOD :          | 390169<br>(99.7%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | GEMFI<br>600MG,STKF2A00 |
| User ID :        | sj97904       | Batch NO. :          | KFM2500260              |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 13 Min 21 Sec   |
| Total Counter :  | 400070        | Total Running Time : | 3296 Hour 45 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 190               | %       | Manual ;           | 45           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 600000            | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 50                | Min     | Supply Time ;      | 0.6          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 1.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 11.31             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 10.87             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 8.87              | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 8.70              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 8.52              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 6.52              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 6.09              | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |     |            |    |
|----------------------|------------|----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.05       | kN | Min : | 1.90       | kN  | Punch NO.; | 17 |
| SD :                 | 4.80       | %  | Max : | 2.48       | kN  | Punch NO.; | 22 |
| Mean Main Pressure : | 8.66       | kN | Min : | 8.18       | kN  | Punch NO.; | 19 |
| SD :                 | 2.50       | %  | Max : | 9.17       | kN  | Punch NO.; | 51 |
| Filling Depth :      | Increase ; | 75 | Times | Decrease ; | 100 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 28       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.90    | mm  | Pre Thickness ;  | 3.91     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.60    | mm  | Main Thickness ; | 3.61     | mm           |
| Filling Depth :               | 10.42   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.775   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 18      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 19      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 44 |
| Lower Punch :    | Average ; | 0.29 | kN |             |    |
|                  | Max ;     | 1.06 | kN | Max Punch ; | 49 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 400070            | Tablets |       |            |         |
| A.W.C Counter : | 391290            | Tablets |       |            |         |
| HEP :           | 611 (0.2%)        | Tablets | LEP : | 510 (0.1%) | Tablets |
| GOOD :          | 390169<br>(99.7%) | Tablets |       |            |         |