

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO200MGSTVI1A00   |
| User ID :        | pb97604       | Batch NO. :          | VIM2401820            |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 48 Min 31 Sec |
| Total Counter :  | 1506688       | Total Running Time : | 332 Hour 48 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 76        | %       | Manual ;           | 45           | RPM       |
| Filling Cam :                      |        | Cam B     |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |        | 4.5-10mm  |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |        | 1740000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 60        | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 40        | Min     | Supply Time ;      | 0.2          | Sec       |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 13.42     | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 12.91     | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 10.53     | kN      | Yes 10             | Times        |           |
| Ref :                              |        | 10.33     | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 10.12     | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 % | 7.74      | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 7.23      | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 30.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None   |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 0.66       | kN | Min : | 0.64       | kN | Punch NO.; | 21 |
| SD :                 | 1.40       | %  | Max : | 0.69       | kN | Punch NO.; | 15 |
| Mean Main Pressure : | 10.32      | kN | Min : | 9.86       | kN | Punch NO.; | 21 |
| SD :                 | 2.10       | %  | Max : | 10.90      | kN | Punch NO.; | 50 |
| Filling Depth :      | Increase ; | 15 | Times | Decrease ; | 21 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.19    | mm  | Pre Thickness ;  | 4.69     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.90    | mm  | Main Thickness ; | 1.40     | mm           |
| Filling Depth :               | 5.73    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.816   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 38      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 30      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 51 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 65 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 49 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1506688            | Tablets |       |            |         |
| A.W.C Counter : | 1486108            | Tablets |       |            |         |
| HEP :           | 513 (0.0%)         | Tablets | LEP : | 562 (0.0%) | Tablets |
| GOOD :          | 1485033<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO200MGSTVI1A00   |
| User ID :        | pb97604       | Batch NO. :          | VIM2401820            |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 48 Min 31 Sec |
| Total Counter :  | 1506688       | Total Running Time : | 332 Hour 48 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 75        | %       | Manual ;           | 59           | RPM       |
| Filling Cam :                      |        | Cam B     |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |        | 4.5-10mm  |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |        | 1740000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 60        | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 40        | Min     | Supply Time ;      | 0.2          | Sec       |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 14.05     | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 13.51     | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 11.03     | kN      | Yes 10             | Times        |           |
| Ref :                              |        | 10.81     | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 10.59     | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 % | 8.11      | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 7.57      | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 30.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 0.76       | kN | Min : | 0.73       | kN | Punch NO.; | 35 |
| SD :                 | 1.90       | %  | Max : | 0.80       | kN | Punch NO.; | 25 |
| Mean Main Pressure : | 10.90      | kN | Min : | 10.48      | kN | Punch NO.; | 23 |
| SD :                 | 2.00       | %  | Max : | 11.53      | kN | Punch NO.; | 55 |
| Filling Depth :      | Increase ; | 17 | Times | Decrease ; | 22 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.39    | mm  | Pre Thickness ;  | 4.89     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.93    | mm  | Main Thickness ; | 1.43     | mm           |
| Filling Depth :               | 5.92    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.816   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 38      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 30      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 15 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 65 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 49 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1506688            | Tablets |       |            |         |
| A.W.C Counter : | 1486108            | Tablets |       |            |         |
| HEP :           | 513 (0.0%)         | Tablets | LEP : | 562 (0.0%) | Tablets |
| GOOD :          | 1485033<br>(99.9%) | Tablets |       |            |         |