

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG<br>STV11B00 |
| User ID :        | pb97604       | Batch NO. :          | VIM2500240              |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 21 Min 06 Sec   |
| Total Counter :  | 950458        | Total Running Time : | 377 Hour 30 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 89                | %       | Manual ;           | 45           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 1500000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 0.3          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 60                | Min     | Supply Time ;      | 0.3          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 18.82             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 18.09             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 14.76             | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 14.47             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 14.18             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 10.86             | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 10.13             | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.64  | kN    | Min :      | 2.55  | kN    | Punch NO.; | 25 |
|                      |            | 1.10  | %     | Max :      | 2.74  | kN    | Punch NO.; | 9  |
| Mean Main Pressure : |            | 14.70 | kN    | Min :      | 14.22 | kN    | Punch NO.; | 9  |
|                      | SD :       | 1.50  | %     | Max :      | 15.40 | kN    | Punch NO.; | 43 |
| Filling Depth :      | Increase ; | 16    | Times | Decrease ; | 45    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.50     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.49    | mm  | Main Thickness ; | 1.99     | mm           |
| Filling Depth :               | 7.88    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.811   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 62      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 44      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 19 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 48 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 48 |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 950458             | Tablets |       |           |         |
| A.W.C Counter : | 938890             | Tablets |       |           |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 22 (0.0%) | Tablets |
| GOOD :          | 938868<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG<br>STV11B00 |
| User ID :        | pb97604       | Batch NO. :          | VIM2500240              |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 21 Min 06 Sec   |
| Total Counter :  | 950458        | Total Running Time : | 377 Hour 30 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 89        | %       | Manual ;           | 59           | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |        | 8.5-14mm  |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |        | 1500000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80        | Min     | Supply Time ;      | 0.3          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 60        | Min     | Supply Time ;      | 0.3          | Sec       |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 19.66     | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 18.90     | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 15.43     | kN      | Yes 10             | Times        |           |
| Ref :                              |        | 15.12     | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 14.82     | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 % | 11.34     | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 10.59     | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 30.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.45       | kN | Min : | 2.38       | kN | Punch NO.; | 32 |
| SD :                 | 1.60       | %  | Max : | 2.57       | kN | Punch NO.; | 46 |
| Mean Main Pressure : | 15.31      | kN | Min : | 14.61      | kN | Punch NO.; | 47 |
| SD :                 | 2.20       | %  | Max : | 16.30      | kN | Punch NO.; | 30 |
| Filling Depth :      | Increase ; | 13 | Times | Decrease ; | 29 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18                | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.50              | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.53    | mm  | Main Thickness ; | 2.03              | mm           |
| Filling Depth :               | 7.87    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.811   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 62      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 44      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 19 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 12 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 48 |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 950458             | Tablets |       |           |         |
| A.W.C Counter : | 938890             | Tablets |       |           |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 22 (0.0%) | Tablets |
| GOOD :          | 938868<br>(100.0%) | Tablets |       |           |         |