

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | AMIOD 200MG<br>STAI3A00 |
| User ID :        | ms88737       | Batch NO. :          | AIM2500120              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 55 Min 08 Sec   |
| Total Counter :  | 1100342       | Total Running Time : | 488 Hour 4 Min          |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 124               | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 1500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 21.86             | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 20.18             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 17.15             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 16.81             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 16.48             | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 %      | 13.45             | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 11.77             | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 20.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.77       | kN  | Min : | 3.34       | kN  | Punch NO.; | 32 |
| SD :                 | 6.30       | %   | Max : | 4.38       | kN  | Punch NO.; | 37 |
| Mean Main Pressure : | 17.28      | kN  | Min : | 16.03      | kN  | Punch NO.; | 10 |
| SD :                 | 4.80       | %   | Max : | 19.32      | kN  | Punch NO.; | 52 |
| Filling Depth :      | Increase ; | 548 | Times | Decrease ; | 554 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 17      | RPM | Capacity :       | 108120   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 21       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.50    | mm  | Pre Thickness ;  | 1.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.95    | mm  | Main Thickness ; | 0.96     | mm           |
| Filling Depth :               | 6.37    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.811   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 61      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.28 | kN |             |    |
|                  | Max ;     | 0.37 | kN | Max Punch ; | 49 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 24 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 1  |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1100342             | Tablets |       |            |         |
| A.W.C Counter : | 1078655             | Tablets |       |            |         |
| HEP :           | 183 (0.0%)          | Tablets | LEP : | 148 (0.0%) | Tablets |
| GOOD :          | 1078324<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | AMIOD 200MG<br>STAI3A00 |
| User ID :        | ms88737       | Batch NO. :          | AIM2500120              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 55 Min 08 Sec   |
| Total Counter :  | 1100342       | Total Running Time : | 488 Hour 4 Min          |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 99                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 1500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 1.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 21.69             | kN      | No                 |              |           |
| HEP :                              | 20.0 % | 20.03             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 17.02             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 16.69             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 16.35             | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 % | 13.35             | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 11.68             | kN      | No                 |              |           |
| SD limit :                         | 20.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 20.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.95       | kN  | Min : | 1.86       | kN  | Punch NO.; | 47 |
| SD :                 | 3.00       | %   | Max : | 2.10       | kN  | Punch NO.; | 6  |
| Mean Main Pressure : | 16.49      | kN  | Min : | 15.49      | kN  | Punch NO.; | 45 |
| SD :                 | 3.00       | %   | Max : | 18.39      | kN  | Punch NO.; | 52 |
| Filling Depth :      | Increase ; | 342 | Times | Decrease ; | 276 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 17      | RPM | Capacity :       | 108120   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 17       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.00    | mm  | Pre Thickness ;  | 2.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.02    | mm  | Main Thickness ; | 1.03     | mm           |
| Filling Depth :               | 6.29    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.811   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 61      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.28 | kN |             |    |
|                  | Max ;     | 0.37 | kN | Max Punch ; | 49 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 24 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 1  |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1100342             | Tablets |       |            |         |
| A.W.C Counter : | 1078655             | Tablets |       |            |         |
| HEP :           | 183 (0.0%)          | Tablets | LEP : | 148 (0.0%) | Tablets |
| GOOD :          | 1078324<br>(100.0%) | Tablets |       |            |         |