

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO800mg<br>STV11C00 |
| User ID :        | ka98463       | Batch NO. :          | VIM2500280              |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 31 Min 07 Sec   |
| Total Counter :  | 300094        | Total Running Time : | 500 Hour 16 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 94       | %       | Manual ;           | 35           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 750000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 25.43    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 24.45    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 19.95    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 19.56    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 19.17    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 14.67    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 13.69    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |              |       |            |       |       |            |    |
|----------------------|--------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.45         | kN    | Min :      | 2.34  | kN    | Punch NO.; | 23 |
| SD :                 | 2.40         | %     | Max :      | 2.58  | kN    | Punch NO.; | 53 |
| Mean Main Pressure : | 19.49        | kN    | Min :      | 18.95 | kN    | Punch NO.; | 31 |
| SD :                 | 1.50         | %     | Max :      | 20.25 | kN    | Punch NO.; | 13 |
| Filling Depth :      | Increase ; 9 | Times | Decrease ; | 75    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 19       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.40    | mm  | Pre Thickness ;  | 4.41     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.76    | mm  | Main Thickness ; | 2.77     | mm           |
| Filling Depth :               | 10.43   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.800   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 77      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.31 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 20 |
| Ejection Force : | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 27 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 300094            | Tablets |       |            |         |
| A.W.C Counter : | 290867            | Tablets |       |            |         |
| HEP :           | 38 (0.0%)         | Tablets | LEP : | 206 (0.1%) | Tablets |
| GOOD :          | 290623<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO800mg<br>STV11C00 |
| User ID :        | ka98463       | Batch NO. :          | VIM2500280              |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 31 Min 07 Sec   |
| Total Counter :  | 300094        | Total Running Time : | 500 Hour 16 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 95                | %       | Manual ;           | 20    | RPM                  |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min : | 6.50 mm              |
| Target Quantity :                  |        | 750000            | Tablets |                    | Max : | 12.00 mm             |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 1.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 25.07             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 24.10             | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 19.67             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 19.28             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 18.90             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 14.46             | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 13.50             | kN      | No                 |       |                      |
| SD limit :                         | 20.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 20.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.91       | kN | Min : | 2.80       | kN | Punch NO.; | 18 |
| SD :                 | 1.80       | %  | Max : | 3.02       | kN | Punch NO.; | 49 |
| Mean Main Pressure : | 19.28      | kN | Min : | 18.29      | kN | Punch NO.; | 28 |
| SD :                 | 1.90       | %  | Max : | 20.08      | kN | Punch NO.; | 27 |
| Filling Depth :      | Increase ; | 11 | Times | Decrease ; | 14 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 19       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.40    | mm  | Pre Thickness ;  | 4.41     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.82    | mm  | Main Thickness ; | 2.83     | mm           |
| Filling Depth :               | 10.50   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.800   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 77      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.31 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 20 |
| Ejection Force : | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 27 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 300094            | Tablets |       |            |         |
| A.W.C Counter : | 290867            | Tablets |       |            |         |
| HEP :           | 38 (0.0%)         | Tablets | LEP : | 206 (0.1%) | Tablets |
| GOOD :          | 290623<br>(99.9%) | Tablets |       |            |         |