

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Mebeve135mg<br>STMV2A00 |
| User ID :        | yu98249       | Batch NO. :          | MVM2500313              |
| Print Interval : | 50000 Tablets | Running Time :       | 10 Hour 14 Min 31 Sec   |
| Total Counter :  | 1300558       | Total Running Time : | 570 Hour 53 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 85       | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 13.81    | kN      | No                 |              |           |
| HEP :                              | 22.0 %      | 12.96    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 10.83    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 10.62    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 10.41    | kN      | Yes 15             | Times        |           |
| LEP :                              | 22.0 %      | 8.28     | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 7.44     | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |      |       |            |      |            |    |
|----------------------|------------|------|-------|------------|------|------------|----|
| Mean Pre Pressure :  | 1.54       | kN   | Min : | 1.48       | kN   | Punch NO.; | 27 |
| SD :                 | 1.70       | %    | Max : | 1.63       | kN   | Punch NO.; | 44 |
| Mean Main Pressure : | 10.59      | kN   | Min : | 10.21      | kN   | Punch NO.; | 22 |
| SD :                 | 1.90       | %    | Max : | 11.12      | kN   | Punch NO.; | 20 |
| Filling Depth :      | Increase ; | 1257 | Times | Decrease ; | 1090 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20     | RPM | Capacity :       | 127200            | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 17                | RPM          |
| Pre Upper Punch Position ;    | 2.50   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.65   | mm  | Pre Thickness ;  | 3.15              | mm           |
| Main Upper Punch Position ;   | 2.20   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.03   | mm  | Main Thickness ; | 1.83              | mm           |
| Filling Depth :               | 9.25   | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.813  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 1      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 8      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.28 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 28 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1300558            | Tablets |       |            |         |
| A.W.C Counter : | 1292064            | Tablets |       |            |         |
| HEP :           | 191 (0.0%)         | Tablets | LEP : | 672 (0.1%) | Tablets |
| GOOD :          | 1291201<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Mebeve135mg<br>STMV2A00 |
| User ID :        | yu98249       | Batch NO. :          | MVM2500313              |
| Print Interval : | 50000 Tablets | Running Time :       | 10 Hour 14 Min 31 Sec   |
| Total Counter :  | 1300558       | Total Running Time : | 570 Hour 53 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 85        | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |        | 8.5-14mm  |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |        | 2000000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 30        | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20        | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 1.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 17.30     | kN      | No                 |              |           |
| HEP :                              | 22.0 % | 16.24     | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 13.58     | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 13.31     | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 13.04     | kN      | Yes 15             | Times        |           |
| LEP :                              | 22.0 % | 10.38     | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 9.32      | kN      | No                 |              |           |
| SD limit :                         | 20.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 20.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.48       | kN  | Min : | 1.41       | kN  | Punch NO.; | 13 |
| SD :                 | 2.20       | %   | Max : | 1.61       | kN  | Punch NO.; | 44 |
| Mean Main Pressure : | 13.31      | kN  | Min : | 12.53      | kN  | Punch NO.; | 28 |
| SD :                 | 2.10       | %   | Max : | 14.35      | kN  | Punch NO.; | 4  |
| Filling Depth :      | Increase ; | 790 | Times | Decrease ; | 731 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20     | RPM | Capacity :       | 127200            | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 17                | RPM          |
| Pre Upper Punch Position ;    | 2.50   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.65   | mm  | Pre Thickness ;  | 3.15              | mm           |
| Main Upper Punch Position ;   | 2.20   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.08   | mm  | Main Thickness ; | 1.88              | mm           |
| Filling Depth :               | 8.51   | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.813  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 1      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 8      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.28 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 28 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1300558            | Tablets |       |            |         |
| A.W.C Counter : | 1292064            | Tablets |       |            |         |
| HEP :           | 191 (0.0%)         | Tablets | LEP : | 672 (0.1%) | Tablets |
| GOOD :          | 1291201<br>(99.9%) | Tablets |       |            |         |