

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Pentoxifylline 400mg
User ID :	pb97604	Batch NO. :	PEN TAB XV 25 001
Print Interval :	50000 Tablets	Running Time :	02 Hour 07 Min 41 Sec
Total Counter :	251438	Total Running Time :	578 Hour 6 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	161	%	Manual ;	15	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		600000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	24.27	kN	No		
HEP :	25.0 %	23.33	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	19.04	kN	Yes 15	Times	
Ref :		18.67	kN			
LCP :	2.0 %	18.29	kN	Yes 15	Times	
LEP :	25.0 %	14.00	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	13.07	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	2.74	kN	Min :	2.36	kN	Punch NO.;	12
		7.70	%	Max :	3.55	kN	Punch NO.;	52
Mean Main Pressure :		18.33	kN	Min :	16.18	kN	Punch NO.;	12
	SD :	5.90	%	Max :	22.13	kN	Punch NO.;	52
Filling Depth :	Increase ;	57	Times	Decrease ;	67	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	127200	Tablets/hour
Feeder : Status ;	Auto		Speed ;	32	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	4.80	mm	Pre Thickness ;	1.81	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.04	mm	Main Thickness ;	1.05	mm
Filling Depth :	8.17	mm	Current Cam :	Cam C	
Main Air pressure :	0.810	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	30	Min			
Lower Punch Remain Time (S2):	3	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.09	kN		
	Max ;	0.21	kN	Max Punch ;	36
Lower Punch :	Average ;	0.20	kN		
	Max ;	0.32	kN	Max Punch ;	15
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	251438	Tablets			
A.W.C Counter :	238793	Tablets			
HEP :	68 (0.0%)	Tablets	LEP :	32 (0.0%)	Tablets
GOOD :	238693 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Pentoxifylline 400mg
User ID :	pb97604	Batch NO. :	PEN TAB XV 25 001
Print Interval :	50000 Tablets	Running Time :	02 Hour 07 Min 41 Sec
Total Counter :	251438	Total Running Time :	578 Hour 6 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	161	%	Manual ;	15	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
Target Quantity :		8.5-14mm			Max :	12.00 mm
Air Pressure Low Limit :		600000	Tablets			
Oil Lubrication :		0.400	MPa			
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	24.95	kN	No		
HEP :	25.0 %	23.99	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	19.57	kN	Yes 15	Times	
Ref :		19.19	kN			
LCP :	2.0 %	18.81	kN	Yes 15	Times	
LEP :	25.0 %	14.39	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	13.43	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	2.65	kN	Min :	2.15	kN	Punch NO.;	23
		7.10	%	Max :	3.05	kN	Punch NO.;	8
Mean Main Pressure :		19.25	kN	Min :	16.43	kN	Punch NO.;	23
	SD :	5.30	%	Max :	21.16	kN	Punch NO.;	8
Filling Depth :	Increase ;	57	Times	Decrease ;	62	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	127200	Tablets/hour
Feeder : Status ;	Auto		Speed ;	32	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	4.80	mm	Pre Thickness ;	1.81	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	4.12	mm	Main Thickness ;	1.12	mm
Filling Depth :	8.12	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.810	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	30	Min			
Lower Punch Remain Time (S2):	3	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.09	kN		
	Max ;	0.22	kN	Max Punch ;	35
Lower Punch :	Average ;	0.20	kN		
	Max ;	0.32	kN	Max Punch ;	15
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	251438	Tablets			
A.W.C Counter :	238793	Tablets			
HEP :	68 (0.0%)	Tablets	LEP :	32 (0.0%)	Tablets
GOOD :	238693 (100.0%)	Tablets			