

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG,STV11B0   |
| User ID :        | pb97604       | Batch NO. :          | VIM2500640            |
| Print Interval : | 50000 Tablets | Running Time :       | 06 Hour 06 Min 30 Sec |
| Total Counter :  | 950366        | Total Running Time : | 596 Hour 51 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |       |                      |
|------------------------------------|--------|-----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 90        | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min : | 6.50 mm              |
|                                    |        | 8.5-14mm  |         |                    | Max : | 12.00 mm             |
| Target Quantity :                  |        | 1500000   | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |           |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30        | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 15        | Min     | Supply Time ;      | 10.0  | Sec                  |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2         | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 1.00      | kN      |                    |       |                      |
|                                    |        |           |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 17.15     | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 16.49     | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 13.45     | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 13.19     | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 12.93     | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 9.89      | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 9.23      | kN      | No                 |       |                      |
| SD limit :                         | 20.0 % |           |         | Yes                |       |                      |
| Pre HSP                            |        | 20.00     | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |           |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.21       | kN | Min : | 2.14       | kN | Punch NO.; | 18 |
| SD :                 | 1.50       | %  | Max : | 2.30       | kN | Punch NO.; | 41 |
| Mean Main Pressure : | 13.10      | kN | Min : | 12.72      | kN | Punch NO.; | 18 |
| SD :                 | 1.10       | %  | Max : | 13.50      | kN | Punch NO.; | 41 |
| Filling Depth :      | Increase ; | 6  | Times | Decrease ; | 0  | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 3.31              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.13    | mm  | Main Thickness ; | 2.14              | mm           |
| Filling Depth :               | 7.14    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.789   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 5       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.25 | kN | Max Punch ; | 27 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 26 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 950366             | Tablets |       |        |         |
| A.W.C Counter : | 949490             | Tablets |       |        |         |
| HEP :           | 1 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 949489<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG,STV11B0   |
| User ID :        | pb97604       | Batch NO. :          | VIM2500640            |
| Print Interval : | 50000 Tablets | Running Time :       | 06 Hour 06 Min 30 Sec |
| Total Counter :  | 950366        | Total Running Time : | 596 Hour 51 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 90        | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |        | 8.5-14mm  |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |        | 1500000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 30        | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 15        | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 1.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 19.37     | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 18.62     | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 15.20     | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 14.90     | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 14.60     | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 11.17     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 10.43     | kN      | No                 |              |           |
| SD limit :                         | 20.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 20.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |              |       |            |       |       |            |    |
|----------------------|--------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.07         | kN    | Min :      | 1.97  | kN    | Punch NO.; | 38 |
| SD :                 | 1.90         | %     | Max :      | 2.17  | kN    | Punch NO.; | 50 |
| Mean Main Pressure : | 15.02        | kN    | Min :      | 14.43 | kN    | Punch NO.; | 38 |
| SD :                 | 1.40         | %     | Max :      | 15.55 | kN    | Punch NO.; | 42 |
| Filling Depth :      | Increase ; 2 | Times | Decrease ; | 0     | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 3.31              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.23    | mm  | Main Thickness ; | 2.24              | mm           |
| Filling Depth :               | 7.37    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.789   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 5       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.25 | kN | Max Punch ; | 27 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 26 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 950366             | Tablets |       |        |         |
| A.W.C Counter : | 949490             | Tablets |       |        |         |
| HEP :           | 1 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 949489<br>(100.0%) | Tablets |       |        |         |