

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SETRA 50<br>MG,STSH4A00 |
| User ID :        | sj97904       | Batch NO. :          | SHM2500350              |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 27 Min 20 Sec   |
| Total Counter :  | 1150878       | Total Running Time : | 614 Hour 24 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 21                | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 6000000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15                | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 17.37             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 16.70             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 13.63             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 13.36             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 13.10             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 10.02             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 9.35              | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 20.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |               |       |            |       |       |            |    |
|----------------------|---------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 0.90          | kN    | Min :      | 0.86  | kN    | Punch NO.; | 27 |
| SD :                 | 1.80          | %     | Max :      | 0.95  | kN    | Punch NO.; | 50 |
| Mean Main Pressure : | 13.31         | kN    | Min :      | 11.99 | kN    | Punch NO.; | 27 |
| SD :                 | 3.90          | %     | Max :      | 14.58 | kN    | Punch NO.; | 50 |
| Filling Depth :      | Increase ; 84 | Times | Decrease ; | 137   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 6        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.80    | mm  | Pre Thickness ;  | 2.80     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.44    | mm  | Main Thickness ; | 1.44     | mm           |
| Filling Depth :               | 5.10    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.769   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 30      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 4       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 25 |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 62 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 13 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1150878             | Tablets |       |        |         |
| A.W.C Counter : | 1147457             | Tablets |       |        |         |
| HEP :           | 2 (0.0%)            | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1147455<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SETRA 50<br>MG,STSH4A00 |
| User ID :        | sj97904       | Batch NO. :          | SHM2500350              |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 27 Min 20 Sec   |
| Total Counter :  | 1150878       | Total Running Time : | 614 Hour 24 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 23                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 6000000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 15                | Min     | Supply Time ;      | 10.0  | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 1.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 18.50             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 17.79             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 14.52             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 14.23             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 13.95             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 10.67             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 9.96              | kN      | No                 |       |                      |
| SD limit :                         | 20.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 20.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.01       | kN | Min : | 0.97       | kN | Punch NO.; | 27 |
| SD :                 | 1.90       | %  | Max : | 1.06       | kN | Punch NO.; | 76 |
| Mean Main Pressure : | 14.12      | kN | Min : | 13.00      | kN | Punch NO.; | 15 |
| SD :                 | 3.90       | %  | Max : | 15.94      | kN | Punch NO.; | 1  |
| Filling Depth :      | Increase ; | 55 | Times | Decrease ; | 70 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 6        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.70    | mm  | Pre Thickness ;  | 2.70     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.57    | mm  | Main Thickness ; | 1.57     | mm           |
| Filling Depth :               | 5.40    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.769   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 30      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 4       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 25 |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 62 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 13 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1150878             | Tablets |       |        |         |
| A.W.C Counter : | 1147457             | Tablets |       |        |         |
| HEP :           | 2 (0.0%)            | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1147455<br>(100.0%) | Tablets |       |        |         |