

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SETRA 50<br>MG,STSH4A00 |
| User ID :        | yu98249       | Batch NO. :          | SHM2500350              |
| Print Interval : | 50000 Tablets | Running Time :       | 17 Hour 03 Min 47 Sec   |
| Total Counter :  | 4418242       | Total Running Time : | 627 Hour 1 Min          |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |       |                      |
|------------------------------------|-------------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ;      | 22                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |             | 6000000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |             |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |             | 30                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |             | 15                | Min     | Supply Time ;      | 10.0  | Sec                  |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |             | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |             | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |             | 1.00              | kN      |                    |       |                      |
|                                    |             |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 %      | 18.09             | kN      | No                 |       |                      |
| HEP :                              | 25.0 %      | 17.39             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %       | 14.19             | kN      | Yes                | 15    | Times                |
| Ref :                              |             | 13.91             | kN      |                    |       |                      |
| LCP :                              | 2.0 %       | 13.64             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 %      | 10.44             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 %      | 9.74              | kN      | No                 |       |                      |
| SD limit :                         | 20.0 %      |                   |         | Yes                |       |                      |
| Pre HSP                            |             | 20.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None        |                   |         |                    |       |                      |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.07  | kN    | Min :      | 1.02  | kN    | Punch NO.; | 22 |
|                      |            | 2.00  | %     | Max :      | 1.15  | kN    | Punch NO.; | 35 |
| Mean Main Pressure : |            | 14.05 | kN    | Min :      | 13.02 | kN    | Punch NO.; | 73 |
|                      | SD :       | 4.50  | %     | Max :      | 16.49 | kN    | Punch NO.; | 35 |
| Filling Depth :      | Increase ; | 316   | Times | Decrease ; | 351   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 6        | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.80   | mm  | Pre Thickness ;  | 2.81     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.49   | mm  | Main Thickness ; | 1.50     | mm           |
| Filling Depth :               | 5.06   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.740  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 17     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 14     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 29 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 62 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 16 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 4418242             | Tablets |       |           |         |
| A.W.C Counter : | 4413633             | Tablets |       |           |         |
| HEP :           | 168 (0.0%)          | Tablets | LEP : | 99 (0.0%) | Tablets |
| GOOD :          | 4413366<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SETRA 50<br>MG,STSH4A00 |
| User ID :        | yu98249       | Batch NO. :          | SHM2500350              |
| Print Interval : | 50000 Tablets | Running Time :       | 17 Hour 03 Min 47 Sec   |
| Total Counter :  | 4418242       | Total Running Time : | 627 Hour 1 Min          |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 25                | %       | Manual ;           | 15           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 6000000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 30                | Min     | Supply Time ;      | 0.5          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 15                | Min     | Supply Time ;      | 10.0         | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 1.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 19.06             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 18.33             | kN      | Yes                | 10 Rotations | 15 Tablets |
| HCP :                              | 2.0 %       | 14.96             | kN      | Yes                | 15 Times     |            |
| Ref :                              |             | 14.67             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 14.37             | kN      | Yes                | 15 Times     |            |
| LEP :                              | 25.0 %      | 11.00             | kN      | Yes                | 10 Rotations | 15 Tablets |
| LSP :                              | 30.0 %      | 10.27             | kN      | No                 |              |            |
| SD limit :                         | 20.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 20.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.18       | kN  | Min : | 1.13       | kN  | Punch NO.; | 19 |
| SD :                 | 2.30       | %   | Max : | 1.30       | kN  | Punch NO.; | 76 |
| Mean Main Pressure : | 14.71      | kN  | Min : | 13.47      | kN  | Punch NO.; | 5  |
| SD :                 | 4.70       | %   | Max : | 17.42      | kN  | Punch NO.; | 76 |
| Filling Depth :      | Increase ; | 138 | Times | Decrease ; | 176 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 7        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.70    | mm  | Pre Thickness ;  | 2.71     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.62    | mm  | Main Thickness ; | 1.62     | mm           |
| Filling Depth :               | 5.47    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.740   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 17      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 14      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 29 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 62 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 16 |

### Tablet Data (N/A)

|                 |                  |         |       |           |         |
|-----------------|------------------|---------|-------|-----------|---------|
| Total Counter : | 4418242          | Tablets |       |           |         |
| A.W.C Counter : | 4413633          | Tablets |       |           |         |
| HEP :           | 168 (0.0%)       | Tablets | LEP : | 99 (0.0%) | Tablets |
| GOOD :          | 4413366 (100.0%) | Tablets |       |           |         |