

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET1000mg STYQ1A00    |
| User ID :        | yu98249       | Batch NO. :          | YQM2501210            |
| Print Interval : | 50000 Tablets | Running Time :       | 14 Hour 11 Min 50 Sec |
| Total Counter :  | 1800706       | Total Running Time : | 706 Hour 5 Min        |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 115      | %       | Manual ;           | 40           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2250000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 27.78    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 26.71    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 21.80    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 21.37    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 20.94    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 16.03    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 14.96    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.27  | kN    | Min :      | 3.03  | kN    | Punch NO.; | 31 |
|                      |            | 3.00  | %     | Max :      | 3.55  | kN    | Punch NO.; | 41 |
| Mean Main Pressure : |            | 21.24 | kN    | Min :      | 20.18 | kN    | Punch NO.; | 4  |
|                      | SD :       | 2.50  | %     | Max :      | 22.57 | kN    | Punch NO.; | 41 |
| Filling Depth :      | Increase ; | 250   | Times | Decrease ; | 264   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20     | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 23       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.70   | mm  | Pre Thickness ;  | 4.71     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.64   | mm  | Main Thickness ; | 3.65     | mm           |
| Filling Depth :               | 9.09   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.747  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 42     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 12     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 17 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 13 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1800706             | Tablets |       |            |         |
| A.W.C Counter : | 1792806             | Tablets |       |            |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 113 (0.0%) | Tablets |
| GOOD :          | 1792693<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET1000mg STYQ1A00    |
| User ID :        | yu98249       | Batch NO. :          | YQM2501210            |
| Print Interval : | 50000 Tablets | Running Time :       | 14 Hour 11 Min 50 Sec |
| Total Counter :  | 1800706       | Total Running Time : | 706 Hour 5 Min        |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 123      | %       | Manual ;           | 40           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2250000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 27.95    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 26.87    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 21.93    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 21.50    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 21.07    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 16.12    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 15.05    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.30  | kN    | Min :      | 3.12  | kN    | Punch NO.; | 46 |
|                      |            | 3.20  | %     | Max :      | 3.51  | kN    | Punch NO.; | 19 |
| Mean Main Pressure : |            | 21.38 | kN    | Min :      | 20.43 | kN    | Punch NO.; | 52 |
|                      | SD :       | 2.50  | %     | Max :      | 22.49 | kN    | Punch NO.; | 8  |
| Filling Depth :      | Increase ; | 118   | Times | Decrease ; | 130   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20     | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 25       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.70   | mm  | Pre Thickness ;  | 4.71     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.72   | mm  | Main Thickness ; | 3.73     | mm           |
| Filling Depth :               | 9.53   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.747  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 42     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 12     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 17 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 13 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1800706             | Tablets |       |            |         |
| A.W.C Counter : | 1792859             | Tablets |       |            |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 113 (0.0%) | Tablets |
| GOOD :          | 1792746<br>(100.0%) | Tablets |       |            |         |