

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO400MG STIU2B00
User ID :	pb97604	Batch NO. :	IUM2500140
Print Interval :	50000 Tablets	Running Time :	00 Hour 29 Min 21 Sec
Total Counter :	50018	Total Running Time :	710 Hour 23 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	200	%	Manual ;	66	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2025000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		60	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	0.5	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	21.68	kN	No		
HEP :	25.0 %	20.85	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	17.01	kN	Yes 15	Times	
Ref :		16.68	kN			
LCP :	2.0 %	16.35	kN	Yes 15	Times	
LEP :	25.0 %	12.51	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	11.68	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	0.51	kN	Min :	0.48	kN	Punch NO.;	6
SD :	2.40	%	Max :	0.53	kN	Punch NO.;	29
Mean Main Pressure :	17.19	kN	Min :	16.00	kN	Punch NO.;	24
SD :	3.40	%	Max :	19.00	kN	Punch NO.;	41
Filling Depth :	Increase ;	336	Times	Decrease ;	327	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	50	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	6.90	mm	Pre Thickness ;	3.91	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.55	mm	Main Thickness ;	1.56	mm
Filling Depth :	7.22	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.800	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	41	Min			
Lower Punch Remain Time (S2):	11	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.11	kN		
	Max ;	0.16	kN	Max Punch ;	3
Ejection Force :	Average ;	0.03	kN		
	Max ;	0.03	kN	Max Punch ;	15

Tablet Data (N/A)

Total Counter :	50018	Tablets			
A.W.C Counter :	37443	Tablets			
HEP :	1 (0.0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	37442 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO400MG STIU2B00
User ID :	pb97604	Batch NO. :	IUM2500140
Print Interval :	50000 Tablets	Running Time :	00 Hour 29 Min 21 Sec
Total Counter :	50018	Total Running Time :	710 Hour 23 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	200	%	Manual ;	74	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2025000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		60	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	0.5	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	21.53	kN	No		
HEP :	25.0 %	20.70	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	16.89	kN	Yes 15	Times	
Ref :		16.56	kN			
LCP :	2.0 %	16.23	kN	Yes 15	Times	
LEP :	25.0 %	12.42	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	11.59	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	0.67	kN	Min :	0.64	kN	Punch NO.;	45
SD :	2.20	%	Max :	0.72	kN	Punch NO.;	20
Mean Main Pressure :	16.86	kN	Min :	15.67	kN	Punch NO.;	52
SD :	4.00	%	Max :	18.86	kN	Punch NO.;	20
Filling Depth :	Increase ;	164	Times	Decrease ;	166	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	50	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	6.90	mm	Pre Thickness ;	3.91	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.63	mm	Main Thickness ;	1.64	mm
Filling Depth :	7.12	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.800	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	41	Min			
Lower Punch Remain Time (S2):	11	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.01	kN	Max Punch ;	32
Lower Punch :	Average ;	0.11	kN		
	Max ;	0.16	kN	Max Punch ;	3
Ejection Force :	Average ;	0.03	kN		
	Max ;	0.03	kN	Max Punch ;	15

Tablet Data (N/A)

Total Counter :	50018	Tablets			
A.W.C Counter :	37443	Tablets			
HEP :	1 (0.0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	37442 (100.0%)	Tablets			