

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Amioda 200mgSTAl3A00
User ID :	am98755	Batch NO. :	AIM2500220
Print Interval :	50000 Tablets	Running Time :	07 Hour 37 Min 06 Sec
Total Counter :	851848	Total Running Time :	731 Hour 48 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	96	%	Manual ;	20	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
Target Quantity :		4.5-10mm			Max :	8.00 mm
Air Pressure Low Limit :		1500000	Tablets			
Oil Lubrication :		0.400	MPa			
Upper Punch Interval Time(S1) ;		25	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	0.5	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	25.60	kN	No		
HEP :	25.0 %	24.62	kN	Yes	10	Rotations 15 Tablets
HCP :	2.5 %	20.19	kN	Yes	15	Times
Ref :		19.70	kN			
LCP :	2.5 %	19.20	kN	Yes	15	Times
LEP :	25.0 %	14.77	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	13.79	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.86	kN	Min :	2.62	kN	Punch NO.;	33
SD :	3.70	%	Max :	3.12	kN	Punch NO.;	45
Mean Main Pressure :	19.44	kN	Min :	18.31	kN	Punch NO.;	18
SD :	2.60	%	Max :	20.80	kN	Punch NO.;	45
Filling Depth :	Increase ;	327	Times	Decrease ;	343	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	17	RPM	Capacity :	108120	Tablets/hour
Feeder : Status ;	Auto		Speed ;	16	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	4.69	mm	Pre Thickness ;	1.70	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	3.83	mm	Main Thickness ;	0.84	mm
Filling Depth :	5.70	mm	Current Cam :	Cam B	
Main Air pressure :	0.795	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	23	Min			
Lower Punch Remain Time (S2):	13	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.02	kN		
	Max ;	0.07	kN	Max Punch ;	1
Lower Punch :	Average ;	0.09	kN		
	Max ;	0.12	kN	Max Punch ;	49
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.06	kN	Max Punch ;	17

Tablet Data (N/A)

Total Counter :	851848	Tablets			
A.W.C Counter :	834251	Tablets			
HEP :	715 (0.1%)	Tablets	LEP :	558 (0.1%)	Tablets
GOOD :	832978 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Amioda 200mgSTAl3A00
User ID :	am98755	Batch NO. :	AIM2500220
Print Interval :	50000 Tablets	Running Time :	07 Hour 37 Min 06 Sec
Total Counter :	851848	Total Running Time :	731 Hour 48 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	90	%	Manual ;	32	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
		4.5-10mm			Max :	8.00 mm
Target Quantity :		1500000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		25	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	0.5	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	25.47	kN	No		
HEP :	25.0 %	24.49	kN	Yes 10	Rotations 15	Tablets
HCP :	2.5 %	20.08	kN	Yes 15	Times	
Ref :		19.59	kN			
LCP :	2.5 %	19.10	kN	Yes 15	Times	
LEP :	25.0 %	14.70	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	13.72	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	2.77	kN	Min :	2.44	kN	Punch NO.;	30
		5.70	%	Max :	3.21	kN	Punch NO.;	17
Mean Main Pressure :		19.34	kN	Min :	17.80	kN	Punch NO.;	20
	SD :	3.50	%	Max :	20.93	kN	Punch NO.;	49
Filling Depth :	Increase ;	288	Times	Decrease ;	285	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	17	RPM	Capacity :	108120	Tablets/hour
Feeder : Status ;	Auto		Speed ;	15	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	4.70	mm	Pre Thickness ;	1.71	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	3.90	mm	Main Thickness ;	0.91	mm
Filling Depth :	6.02	mm	Current Cam :	Cam B	
Main Air pressure :	0.795	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	23	Min			
Lower Punch Remain Time (S2):	13	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.02	kN		
	Max ;	0.08	kN	Max Punch ;	1
Lower Punch :	Average ;	0.09	kN		
	Max ;	0.12	kN	Max Punch ;	49
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.06	kN	Max Punch ;	17

Tablet Data (N/A)

Total Counter :	851848	Tablets			
A.W.C Counter :	834251	Tablets			
HEP :	715 (0.1%)	Tablets	LEP :	558 (0.1%)	Tablets
GOOD :	832978 (99.8%)	Tablets			