

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Amioda 200mgSTAl3A00  |
| User ID :        | pb97604       | Batch NO. :          | AIM2500320            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 31 Min 31 Sec |
| Total Counter :  | 50004         | Total Running Time : | 737 Hour 51 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 92       | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 19.67    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 18.91    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.5 %       | 15.51    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 15.13    | kN      |                    |              |           |
| LCP :                              | 2.5 %       | 14.75    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 11.35    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 10.59    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.74       | kN  | Min : | 2.50       | kN  | Punch NO.; | 3  |
| SD :                 | 4.10       | %   | Max : | 3.00       | kN  | Punch NO.; | 46 |
| Mean Main Pressure : | 14.95      | kN  | Min : | 13.79      | kN  | Punch NO.; | 3  |
| SD :                 | 3.10       | %   | Max : | 15.97      | kN  | Punch NO.; | 46 |
| Filling Depth :      | Increase ; | 869 | Times | Decrease ; | 828 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 14       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.69    | mm  | Pre Thickness ;  | 1.70     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.00    | mm  | Main Thickness ; | 1.01     | mm           |
| Filling Depth :               | 4.72    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.796   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 24      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 24      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 5  |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 49 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 17 |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 50004             | Tablets |       |        |         |
| A.W.C Counter : | 43396             | Tablets |       |        |         |
| HEP :           | 0 (0%)            | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 43396<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Amioda 200mgSTAI3A00  |
| User ID :        | pb97604       | Batch NO. :          | AIM2500320            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 31 Min 31 Sec |
| Total Counter :  | 50004         | Total Running Time : | 737 Hour 51 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 95       | %       | Manual ;           | 32           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 19.01    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 18.28    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.5 %       | 14.99    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 14.63    | kN      |                    |              |           |
| LCP :                              | 2.5 %       | 14.26    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 10.97    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 10.24    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.90           | kN    | Min :      | 2.57  | kN    | Punch NO.; | 24 |
| SD :                 | 5.40           | %     | Max :      | 3.31  | kN    | Punch NO.; | 17 |
| Mean Main Pressure : | 14.73          | kN    | Min :      | 13.67 | kN    | Punch NO.; | 24 |
| SD :                 | 3.40           | %     | Max :      | 16.07 | kN    | Punch NO.; | 17 |
| Filling Depth :      | Increase ; 597 | Times | Decrease ; | 563   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 14       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.70    | mm  | Pre Thickness ;  | 1.71     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.15    | mm  | Main Thickness ; | 1.16     | mm           |
| Filling Depth :               | 5.34    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.796   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 24      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 24      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 5  |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 49 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 17 |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 50004             | Tablets |       |        |         |
| A.W.C Counter : | 43396             | Tablets |       |        |         |
| HEP :           | 0 (0%)            | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 43396<br>(100.0%) | Tablets |       |        |         |