

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Amioda 200mgSTAI3A00  |
| User ID :        | pn98247       | Batch NO. :          | AIM2500320            |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 10 Min 52 Sec |
| Total Counter :  | 806660        | Total Running Time : | 744 Hour 28 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 92       | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 20.41    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 19.62    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.5 %       | 16.09    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 15.70    | kN      |                    |              |           |
| LCP :                              | 2.5 %       | 15.31    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 11.77    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 10.99    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.94       | kN  | Min : | 2.74       | kN  | Punch NO.; | 48 |
| SD :                 | 4.10       | %   | Max : | 3.23       | kN  | Punch NO.; | 6  |
| Mean Main Pressure : | 15.75      | kN  | Min : | 14.48      | kN  | Punch NO.; | 48 |
| SD :                 | 3.60       | %   | Max : | 17.17      | kN  | Punch NO.; | 6  |
| Filling Depth :      | Increase ; | 257 | Times | Decrease ; | 264 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 17      | RPM | Capacity :       | 108120   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 16       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.69    | mm  | Pre Thickness ;  | 1.70     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.99    | mm  | Main Thickness ; | 1.00     | mm           |
| Filling Depth :               | 5.24    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.812   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 7       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 22      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 4  |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 17 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 28 |

### Tablet Data (N/A)

|                 |                 |         |       |           |         |
|-----------------|-----------------|---------|-------|-----------|---------|
| Total Counter : | 806660          | Tablets |       |           |         |
| A.W.C Counter : | 792649          | Tablets |       |           |         |
| HEP :           | 159 (0.0%)      | Tablets | LEP : | 61 (0.0%) | Tablets |
| GOOD :          | 792429 (100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Amioda 200mgSTAl3A00  |
| User ID :        | pn98247       | Batch NO. :          | AIM2500320            |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 10 Min 52 Sec |
| Total Counter :  | 806660        | Total Running Time : | 744 Hour 28 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 95       | %       | Manual ;           | 32           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 18.93    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 18.20    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.5 %       | 14.93    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 14.56    | kN      |                    |              |           |
| LCP :                              | 2.5 %       | 14.20    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 10.92    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 10.19    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.83       | kN  | Min : | 2.46       | kN  | Punch NO.; | 36 |
| SD :                 | 6.70       | %   | Max : | 3.36       | kN  | Punch NO.; | 6  |
| Mean Main Pressure : | 14.69      | kN  | Min : | 13.61      | kN  | Punch NO.; | 36 |
| SD :                 | 4.10       | %   | Max : | 16.17      | kN  | Punch NO.; | 6  |
| Filling Depth :      | Increase ; | 447 | Times | Decrease ; | 444 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 17      | RPM | Capacity :       | 108120            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 16                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.70    | mm  | Pre Thickness ;  | 1.71              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.14    | mm  | Main Thickness ; | 1.15              | mm           |
| Filling Depth :               | 5.69    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.812   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 7       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 22      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 4  |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 17 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 17 |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 806660             | Tablets |       |           |         |
| A.W.C Counter : | 792649             | Tablets |       |           |         |
| HEP :           | 159 (0.0%)         | Tablets | LEP : | 61 (0.0%) | Tablets |
| GOOD :          | 792429<br>(100.0%) | Tablets |       |           |         |