

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Sertra 100mgSTSH4B00  |
| User ID :        | ms98171       | Batch NO. :          | SHM2500910            |
| Print Interval : | 50000 Tablets | Running Time :       | 11 Hour 02 Min 33 Sec |
| Total Counter :  | 2852270       | Total Running Time : | 811 Hour 19 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 26       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 3000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 40       | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 16.34    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 15.72    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 12.82    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 12.57    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 12.32    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 9.43     | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 8.80     | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |               |       |            |       |       |            |    |
|----------------------|---------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.45          | kN    | Min :      | 2.36  | kN    | Punch NO.; | 17 |
| SD :                 | 1.50          | %     | Max :      | 2.53  | kN    | Punch NO.; | 46 |
| Mean Main Pressure : | 12.49         | kN    | Min :      | 11.71 | kN    | Punch NO.; | 8  |
| SD :                 | 2.90          | %     | Max :      | 13.54 | kN    | Punch NO.; | 75 |
| Filling Depth :      | Increase ; 88 | Times | Decrease ; | 108   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 7        | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.60    | mm  | Pre Thickness ;  | 3.10     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.55    | mm  | Main Thickness ; | 2.05     | mm           |
| Filling Depth :               | 7.58    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.780   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 24      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 24 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.34 | kN | Max Punch ; | 57 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 42 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 2852270             | Tablets |       |        |         |
| A.W.C Counter : | 2843387             | Tablets |       |        |         |
| HEP :           | 10 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 2843377<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                                                                                                                                                                         |                                                                                                                                                                           |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Machine Name :</b> VANTIX P720D<br><b>User ID :</b> ms98171<br><b>Print Interval :</b> 50000 Tablets<br><b>Total Counter :</b> 2852270<br><b>Station :</b> Station 2 | <b>Product Name :</b> Sertra 100mgSTSH4B00<br><b>Batch NO. :</b> SHM2500910<br><b>Running Time :</b> 11 Hour 02 Min 33 Sec<br><b>Total Running Time :</b> 811 Hour 19 Min |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 25       | %       | Manual ;           | 48           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 3000000  | Tablets |                    |              |           |
| Oil Lubrication :                  |             | 0.600    | MPa     |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 40       | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 16.95    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 16.30    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 13.30    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 13.04    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 12.78    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 9.78     | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 9.13     | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.48  | kN    | Min :      | 2.42  | kN    | Punch NO.; | 60 |
|                      |            | 1.60  | %     | Max :      | 2.67  | kN    | Punch NO.; | 75 |
| Mean Main Pressure : |            | 13.09 | kN    | Min :      | 12.37 | kN    | Punch NO.; | 8  |
|                      | SD :       | 3.00  | %     | Max :      | 15.06 | kN    | Punch NO.; | 75 |
| Filling Depth :      | Increase ; | 52    | Times | Decrease ; | 62    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 7        | RPM          |
| Pre Upper Punch Position ;    | 2.50   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.60   | mm  | Pre Thickness ;  | 3.10     | mm           |
| Main Upper Punch Position ;   | 2.50   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.70   | mm  | Main Thickness ; | 2.20     | mm           |
| Filling Depth :               | 7.77   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.780  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 24     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 24 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.34 | kN | Max Punch ; | 57 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 42 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 2852270             | Tablets |       |        |         |
| A.W.C Counter : | 2843464             | Tablets |       |        |         |
| HEP :           | 10 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 2843454<br>(100.0%) | Tablets |       |        |         |