

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | rp98927       | Batch NO. :          | MVM2500413              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 25 Min 34 Sec   |
| Total Counter :  | 51674         | Total Running Time : | 812 Hour 13 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 80                | %       | Manual ;           | 25    | RPM                  |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min : | 6.50 mm              |
| Target Quantity :                  |        | 3000000           | Tablets |                    | Max : | 12.00 mm             |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 40                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 16.34             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 15.72             | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 12.82             | kN      | Yes                | 10    | Times                |
| Ref :                              |        | 12.57             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 12.32             | kN      | Yes                | 10    | Times                |
| LEP :                              | 25.0 % | 9.43              | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 8.80              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.04       | kN | Min : | 2.01       | kN | Punch NO.; | 15 |
| SD :                 | 1.10       | %  | Max : | 2.13       | kN | Punch NO.; | 23 |
| Mean Main Pressure : | 12.50      | kN | Min : | 11.32      | kN | Punch NO.; | 19 |
| SD :                 | 4.30       | %  | Max : | 13.79      | kN | Punch NO.; | 8  |
| Filling Depth :      | Increase ; | 7  | Times | Decrease ; | 7  | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 16       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.50    | mm  | Main Thickness ; | 1.51     | mm           |
| Filling Depth :               | 9.28    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.795   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 40      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 11      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 41 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 23 |
| Ejection Force : | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.36 | kN | Max Punch ; | 23 |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 51674             | Tablets |       |        |         |
| A.W.C Counter : | 47328             | Tablets |       |        |         |
| HEP :           | 3 (0.0%)          | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 47325<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | rp98927       | Batch NO. :          | MVM2500413              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 25 Min 34 Sec   |
| Total Counter :  | 51674         | Total Running Time : | 812 Hour 13 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 80                | %       | Manual ;           | 25           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 3000000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 40                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 20.88             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 20.07             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 16.38             | kN      | Yes 10             | Times        |           |
| Ref :                              |        | 16.06             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 15.74             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 % | 12.04             | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 11.24             | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.83       | kN | Min : | 1.78       | kN | Punch NO.; | 28 |
| SD :                 | 1.40       | %  | Max : | 1.94       | kN | Punch NO.; | 23 |
| Mean Main Pressure : | 15.81      | kN | Min : | 14.03      | kN | Punch NO.; | 11 |
| SD :                 | 4.10       | %  | Max : | 17.50      | kN | Punch NO.; | 44 |
| Filling Depth :      | Increase ; | 13 | Times | Decrease ; | 6  | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 16       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.50    | mm  | Main Thickness ; | 1.50     | mm           |
| Filling Depth :               | 9.20    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.795   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 40      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 11      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.23 | kN | Max Punch ; | 41 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 23 |
| Ejection Force : | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.36 | kN | Max Punch ; | 23 |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 51674             | Tablets |       |        |         |
| A.W.C Counter : | 47328             | Tablets |       |        |         |
| HEP :           | 3 (0.0%)          | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 47325<br>(100.0%) | Tablets |       |        |         |