

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MEBEV 135 mgSTMV2A00
User ID :	pb97604	Batch NO. :	MVM2500413
Print Interval :	50000 Tablets	Running Time :	05 Hour 06 Min 15 Sec
Total Counter :	644838	Total Running Time :	816 Hour 54 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	80	%	Manual ;	25	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		3000000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	17.35	kN	No		
HEP :	25.0 %	16.68	kN	Yes	10	Rotations 10 Tablets
HCP :	2.0 %	13.61	kN	Yes	10	Times
Ref :		13.35	kN			
LCP :	2.0 %	13.08	kN	Yes	10	Times
LEP :	25.0 %	10.01	kN	Yes	10	Rotations 10 Tablets
LSP :	30.0 %	9.34	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.11	kN	Min :	2.07	kN	Punch NO.;	2
SD :	1.30	%	Max :	2.20	kN	Punch NO.;	48
Mean Main Pressure :	13.38	kN	Min :	12.17	kN	Punch NO.;	19
SD :	3.90	%	Max :	14.61	kN	Punch NO.;	44
Filling Depth :	Increase ;	188	Times	Decrease ;	156	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	127200	Tablets/hour
Feeder : Status ;	Auto		Speed ;	16	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.65	mm	Main Thickness ;	1.66	mm
Filling Depth :	9.44	mm	Current Cam :	Cam C	
Main Air pressure :	0.815	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	6	Min			
Lower Punch Remain Time (S2):	7	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.15	kN		
	Max ;	0.22	kN	Max Punch ;	43
Lower Punch :	Average ;	0.11	kN		
	Max ;	0.31	kN	Max Punch ;	48
Ejection Force :	Average ;	0.09	kN		
	Max ;	0.40	kN	Max Punch ;	48

Tablet Data (N/A)

Total Counter :	644838	Tablets			
A.W.C Counter :	639736	Tablets			
HEP :	109 (0.0%)	Tablets	LEP :	13 (0.0%)	Tablets
GOOD :	639614 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MEBEV 135 mgSTMV2A00
User ID :	pb97604	Batch NO. :	MVM2500413
Print Interval :	50000 Tablets	Running Time :	05 Hour 06 Min 15 Sec
Total Counter :	644838	Total Running Time :	816 Hour 54 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	80	%	Manual ;	25	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		3000000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	17.45	kN	No		
HEP :	25.0 %	16.78	kN	Yes 10	Rotations 10	Tablets
HCP :	2.0 %	13.69	kN	Yes 10	Times	
Ref :		13.43	kN			
LCP :	2.0 %	13.16	kN	Yes 10	Times	
LEP :	25.0 %	10.07	kN	Yes 10	Rotations 10	Tablets
LSP :	30.0 %	9.40	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.03	kN	Min :	1.98	kN	Punch NO.;	42
SD :	1.40	%	Max :	2.15	kN	Punch NO.;	48
Mean Main Pressure :	13.26	kN	Min :	12.15	kN	Punch NO.;	37
SD :	3.90	%	Max :	14.34	kN	Punch NO.;	50
Filling Depth :	Increase ;	44	Times	Decrease ;	33	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	127200	Tablets/hour
Feeder : Status ;	Auto		Speed ;	16	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.65	mm	Main Thickness ;	1.66	mm
Filling Depth :	9.27	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.815	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	6	Min			
Lower Punch Remain Time (S2):	7	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.15	kN		
	Max ;	0.23	kN	Max Punch ;	41
Lower Punch :	Average ;	0.11	kN		
	Max ;	0.31	kN	Max Punch ;	48
Ejection Force :	Average ;	0.09	kN		
	Max ;	0.40	kN	Max Punch ;	48

Tablet Data (N/A)

Total Counter :	644838	Tablets			
A.W.C Counter :	639736	Tablets			
HEP :	109 (0.0%)	Tablets	LEP :	13 (0.0%)	Tablets
GOOD :	639614 (100.0%)	Tablets			