

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | pr190041      | Batch NO. :          | MVM2500513              |
| Print Interval : | 50000 Tablets | Running Time :       | 11 Hour 10 Min 42 Sec   |
| Total Counter :  | 1417334       | Total Running Time : | 838 Hour 18 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 102               | %       | Manual ;           | 34           | RPM        |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm    |
| Target Quantity :                  |             | 2000000           | Tablets |                    | Max :        | 12.00 mm   |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 30                | Min     | Supply Time ;      | 0.2          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 16.68             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 16.03             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 13.08             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 12.83             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 12.57             | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 9.62              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 8.98              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | None        |                   |         |                    |              |            |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.24       | kN | Min : | 2.18       | kN | Punch NO.; | 5  |
| SD :                 | 1.80       | %  | Max : | 2.39       | kN | Punch NO.; | 28 |
| Mean Main Pressure : | 12.41      | kN | Min : | 10.94      | kN | Punch NO.; | 53 |
| SD :                 | 4.90       | %  | Max : | 13.71      | kN | Punch NO.; | 2  |
| Filling Depth :      | Increase ; | 82 | Times | Decrease ; | 68 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 20       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.67    | mm  | Main Thickness ; | 1.68     | mm           |
| Filling Depth :               | 9.43    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.770   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 13      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 15      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 43 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.32 | kN | Max Punch ; | 28 |
| Ejection Force : | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.56 | kN | Max Punch ; | 28 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 1417334             | Tablets |       |           |         |
| A.W.C Counter : | 1410630             | Tablets |       |           |         |
| HEP :           | 440 (0.0%)          | Tablets | LEP : | 52 (0.0%) | Tablets |
| GOOD :          | 1410138<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | pr190041      | Batch NO. :          | MVM2500513              |
| Print Interval : | 50000 Tablets | Running Time :       | 11 Hour 10 Min 42 Sec   |
| Total Counter :  | 1417334       | Total Running Time : | 838 Hour 18 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 102       | %       | Manual ;           | 25           | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |        | 8.5-14mm  |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |        | 2000000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.600     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 30        | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20        | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 17.55     | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 16.88     | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 13.77     | kN      | Yes 10             | Times        |           |
| Ref :                              |        | 13.50     | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 13.23     | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 % | 10.13     | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 9.45      | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 30.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.28       | kN | Min : | 2.19       | kN | Punch NO.; | 28 |
| SD :                 | 1.20       | %  | Max : | 2.36       | kN | Punch NO.; | 44 |
| Mean Main Pressure : | 13.72      | kN | Min : | 12.48      | kN | Punch NO.; | 53 |
| SD :                 | 4.90       | %  | Max : | 15.11      | kN | Punch NO.; | 9  |
| Filling Depth :      | Increase ; | 67 | Times | Decrease ; | 65 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 20                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.68    | mm  | Main Thickness ; | 1.69              | mm           |
| Filling Depth :               | 8.97    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.770   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 13      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 15      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 43 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.31 | kN | Max Punch ; | 11 |
| Ejection Force : | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.56 | kN | Max Punch ; | 28 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 1417334             | Tablets |       |           |         |
| A.W.C Counter : | 1410630             | Tablets |       |           |         |
| HEP :           | 440 (0.0%)          | Tablets | LEP : | 52 (0.0%) | Tablets |
| GOOD :          | 1410138<br>(100.0%) | Tablets |       |           |         |