

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | rp98927       | Batch NO. :          | MVM2500613              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 05 Min 30 Sec   |
| Total Counter :  | 11608         | Total Running Time : | 842 Hour 49 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 81                | %       | Manual ;           | 34           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 2000000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.500             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 16.48             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 15.84             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 12.93             | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 12.67             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 12.42             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 9.51              | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 8.87              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.06  | kN    | Min :      | 2.00  | kN    | Punch NO.; | 15 |
|                      |            | 2.80  | %     | Max :      | 2.36  | kN    | Punch NO.; | 36 |
| Mean Main Pressure : |            | 12.75 | kN    | Min :      | 11.61 | kN    | Punch NO.; | 26 |
|                      | SD :       | 3.80  | %     | Max :      | 14.01 | kN    | Punch NO.; | 38 |
| Filling Depth :      | Increase ; | 0     | Times | Decrease ; | 3     | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 16       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.67    | mm  | Main Thickness ; | 1.68     | mm           |
| Filling Depth :               | 8.81    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.780   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 13      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 15      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 17 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.65 | kN | Max Punch ; | 36 |
| Ejection Force : | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 1.03 | kN | Max Punch ; | 36 |

### Tablet Data (N/A)

|                 |                  |         |       |        |         |
|-----------------|------------------|---------|-------|--------|---------|
| Total Counter : | 11608            | Tablets |       |        |         |
| A.W.C Counter : | 9753             | Tablets |       |        |         |
| HEP :           | 0 (0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 9753<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | rp98927       | Batch NO. :          | MVM2500613              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 05 Min 30 Sec   |
| Total Counter :  | 11608         | Total Running Time : | 842 Hour 49 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 81                | %       | Manual ;           | 25           | RPM        |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm    |
| Target Quantity :                  |             | 2000000           | Tablets |                    | Max :        | 12.00 mm   |
| Air Pressure Low Limit :           |             | 0.500             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 30                | Min     | Supply Time ;      | 0.2          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 15.68             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 15.08             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 12.31             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 12.06             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 11.82             | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 9.05              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 8.44              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.89       | kN | Min : | 1.78       | kN | Punch NO.; | 36 |
| SD :                 | 1.20       | %  | Max : | 1.93       | kN | Punch NO.; | 26 |
| Mean Main Pressure : | 12.02      | kN | Min : | 10.07      | kN | Punch NO.; | 49 |
| SD :                 | 5.80       | %  | Max : | 13.78      | kN | Punch NO.; | 23 |
| Filling Depth :      | Increase ; | 1  | Times | Decrease ; | 3  | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 16       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.68    | mm  | Main Thickness ; | 1.69     | mm           |
| Filling Depth :               | 8.83    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.780   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 13      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 15      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 19 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.65 | kN | Max Punch ; | 36 |
| Ejection Force : | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 1.03 | kN | Max Punch ; | 36 |

### Tablet Data (N/A)

|                 |                  |         |       |        |         |
|-----------------|------------------|---------|-------|--------|---------|
| Total Counter : | 11608            | Tablets |       |        |         |
| A.W.C Counter : | 9753             | Tablets |       |        |         |
| HEP :           | 0 (0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 9753<br>(100.0%) | Tablets |       |        |         |