

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | pb97604       | Batch NO. :          | MVM2500813              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 17 Min 50 Sec   |
| Total Counter :  | 1300548       | Total Running Time : | 881 Hour 30 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 52                | %       | Manual ;           | 11           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 2000000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 15.96             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 15.35             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 12.52             | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 12.28             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 12.03             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 9.21              | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 8.59              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.19       | kN  | Min : | 2.15       | kN  | Punch NO.; | 43 |
| SD :                 | 0.70       | %   | Max : | 2.24       | kN  | Punch NO.; | 27 |
| Mean Main Pressure : | 12.26      | kN  | Min : | 11.33      | kN  | Punch NO.; | 19 |
| SD :                 | 3.20       | %   | Max : | 12.96      | kN  | Punch NO.; | 41 |
| Filling Depth :      | Increase ; | 819 | Times | Decrease ; | 752 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22     | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 11       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.67   | mm  | Main Thickness ; | 1.68     | mm           |
| Filling Depth :               | 8.76   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.765  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 1      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 13     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.26 | kN |             |    |
|                  | Max ;     | 0.34 | kN | Max Punch ; | 28 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 19 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                     |         |       |          |         |
|-----------------|---------------------|---------|-------|----------|---------|
| Total Counter : | 1300548             | Tablets |       |          |         |
| A.W.C Counter : | 1297349             | Tablets |       |          |         |
| HEP :           | 25 (0.0%)           | Tablets | LEP : | 6 (0.0%) | Tablets |
| GOOD :          | 1297318<br>(100.0%) | Tablets |       |          |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | pb97604       | Batch NO. :          | MVM2500813              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 17 Min 50 Sec   |
| Total Counter :  | 1300548       | Total Running Time : | 881 Hour 30 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 55                | %       | Manual ;           | 11    | RPM                  |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min : | 6.50 mm              |
| Target Quantity :                  |        | 2000000           | Tablets |                    | Max : | 12.00 mm             |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.2   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 15.83             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 15.22             | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 12.42             | kN      | Yes                | 10    | Times                |
| Ref :                              |        | 12.18             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 11.93             | kN      | Yes                | 10    | Times                |
| LEP :                              | 25.0 % | 9.13              | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 8.52              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.30       | kN  | Min : | 2.26       | kN  | Punch NO.; | 41 |
| SD :                 | 0.70       | %   | Max : | 2.34       | kN  | Punch NO.; | 53 |
| Mean Main Pressure : | 12.28      | kN  | Min : | 11.07      | kN  | Punch NO.; | 53 |
| SD :                 | 5.00       | %   | Max : | 13.43      | kN  | Punch NO.; | 9  |
| Filling Depth :      | Increase ; | 833 | Times | Decrease ; | 843 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 22     | RPM | Capacity :       | 139920            | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.51              | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.68   | mm  | Main Thickness ; | 1.69              | mm           |
| Filling Depth :               | 8.84   | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.765  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 1      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 13     | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.26 | kN |             |    |
|                  | Max ;     | 0.34 | kN | Max Punch ; | 28 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 19 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                     |         |       |          |         |
|-----------------|---------------------|---------|-------|----------|---------|
| Total Counter : | 1300548             | Tablets |       |          |         |
| A.W.C Counter : | 1297349             | Tablets |       |          |         |
| HEP :           | 25 (0.0%)           | Tablets | LEP : | 6 (0.0%) | Tablets |
| GOOD :          | 1297318<br>(100.0%) | Tablets |       |          |         |