

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Prop 40mg STPJ2B00    |
| User ID :        | yu98249       | Batch NO. :          | PJM2500440            |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 02 Min 21 Sec |
| Total Counter :  | 1278980       | Total Running Time : | 919 Hour 7 Min        |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 40       | %       | Manual ;           | 11           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 29.27    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 28.14    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 22.97    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 22.52    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 22.06    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 16.89    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 15.76    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |              |       |              |       |    |            |    |
|----------------------|--------------|-------|--------------|-------|----|------------|----|
| Mean Pre Pressure :  | 1.85         | kN    | Min :        | 1.80  | kN | Punch NO.; | 9  |
| SD :                 | 1.90         | %     | Max :        | 1.93  | kN | Punch NO.; | 12 |
| Mean Main Pressure : | 22.34        | kN    | Min :        | 21.14 | kN | Punch NO.; | 14 |
| SD :                 | 2.50         | %     | Max :        | 23.61 | kN | Punch NO.; | 16 |
| Filling Depth :      | Increase ; 5 | Times | Decrease ; 0 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25     | RPM | Capacity :       | 159000   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.56   | mm  | Main Thickness ; | 1.57     | mm           |
| Filling Depth :               | 5.67   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.784  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 13     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 19     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.35 | kN |             |    |
|                  | Max ;     | 0.41 | kN | Max Punch ; | 39 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 41 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 53 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1278980             | Tablets |       |            |         |
| A.W.C Counter : | 1275864             | Tablets |       |            |         |
| HEP :           | 2 (0.0%)            | Tablets | LEP : | 387 (0.0%) | Tablets |
| GOOD :          | 1275475<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Prop 40mg STPJ2B00    |
| User ID :        | yu98249       | Batch NO. :          | PJM2500440            |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 02 Min 21 Sec |
| Total Counter :  | 1278980       | Total Running Time : | 919 Hour 7 Min        |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 39       | %       | Manual ;           | 11           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 30.05    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 28.89    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 23.58    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 23.11    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 22.65    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 17.34    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 16.18    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |              |       |              |       |    |            |    |
|----------------------|--------------|-------|--------------|-------|----|------------|----|
| Mean Pre Pressure :  | 1.98         | kN    | Min :        | 1.94  | kN | Punch NO.; | 4  |
| SD :                 | 0.80         | %     | Max :        | 2.04  | kN | Punch NO.; | 53 |
| Mean Main Pressure : | 23.35        | kN    | Min :        | 21.79 | kN | Punch NO.; | 1  |
| SD :                 | 3.60         | %     | Max :        | 25.43 | kN | Punch NO.; | 30 |
| Filling Depth :      | Increase ; 3 | Times | Decrease ; 3 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 159000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.55    | mm  | Main Thickness ; | 1.56     | mm           |
| Filling Depth :               | 5.56    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.784   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 13      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 19      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.35 | kN |             |    |
|                  | Max ;     | 0.41 | kN | Max Punch ; | 39 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 41 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 41 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1278980             | Tablets |       |            |         |
| A.W.C Counter : | 1275864             | Tablets |       |            |         |
| HEP :           | 2 (0.0%)            | Tablets | LEP : | 387 (0.0%) | Tablets |
| GOOD :          | 1275475<br>(100.0%) | Tablets |       |            |         |