

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Prop 40mg STPJ2B00    |
| User ID :        | yu98249       | Batch NO. :          | PJM2500640            |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 57 Min 54 Sec |
| Total Counter :  | 1585878       | Total Running Time : | 945 Hour 50 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |       |                      |
|------------------------------------|-------------|----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ;      | 42       | %       | Manual ;           | 11    | RPM                  |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min : | 2.50 mm              |
|                                    |             | 4.5-10mm |         |                    | Max : | 8.00 mm              |
| Target Quantity :                  |             | 2000000  | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |       |                      |
| Oil Lubrication :                  |             |          |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |             | 30       | Min     | Supply Time ;      | 0.2   | Sec                  |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |       |                      |
| Initial Reject Time :              |             | 2        | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |       |                      |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |       |                      |
|                                    |             |          |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 %      | 30.75    | kN      | No                 |       |                      |
| HEP :                              | 25.0 %      | 29.57    | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %       | 24.13    | kN      | Yes                | 10    | Times                |
| Ref :                              |             | 23.66    | kN      |                    |       |                      |
| LCP :                              | 2.0 %       | 23.18    | kN      | Yes                | 10    | Times                |
| LEP :                              | 25.0 %      | 17.74    | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 %      | 16.56    | kN      | No                 |       |                      |
| SD limit :                         | 35.0 %      |          |         | Yes                |       |                      |
| Pre HSP                            |             | 30.00    | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None        |          |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.88       | kN | Min : | 1.83       | kN | Punch NO.; | 32 |
| SD :                 | 1.20       | %  | Max : | 1.93       | kN | Punch NO.; | 49 |
| Mean Main Pressure : | 23.73      | kN | Min : | 22.34      | kN | Punch NO.; | 53 |
| SD :                 | 2.80       | %  | Max : | 25.29      | kN | Punch NO.; | 16 |
| Filling Depth :      | Increase ; | 0  | Times | Decrease ; | 19 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 159000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.56    | mm  | Main Thickness ; | 1.57     | mm           |
| Filling Depth :               | 5.66    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.777   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 10      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 18      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.36 | kN |             |    |
|                  | Max ;     | 0.43 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 41 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 41 |

### Tablet Data (N/A)

|                 |                     |         |       |          |         |
|-----------------|---------------------|---------|-------|----------|---------|
| Total Counter : | 1585878             | Tablets |       |          |         |
| A.W.C Counter : | 1581731             | Tablets |       |          |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 5 (0.0%) | Tablets |
| GOOD :          | 1581726<br>(100.0%) | Tablets |       |          |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Prop 40mg STPJ2B00    |
| User ID :        | yu98249       | Batch NO. :          | PJM2500640            |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 57 Min 54 Sec |
| Total Counter :  | 1585878       | Total Running Time : | 945 Hour 50 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 39       | %       | Manual ;           | 11           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 31.53    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 30.32    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 24.74    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 24.26    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 23.77    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 18.19    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 16.98    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |              |       |              |       |    |            |    |
|----------------------|--------------|-------|--------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.04         | kN    | Min :        | 2.00  | kN | Punch NO.; | 32 |
| SD :                 | 1.10         | %     | Max :        | 2.10  | kN | Punch NO.; | 38 |
| Mean Main Pressure : | 24.10        | kN    | Min :        | 21.90 | kN | Punch NO.; | 53 |
| SD :                 | 3.60         | %     | Max :        | 26.16 | kN | Punch NO.; | 2  |
| Filling Depth :      | Increase ; 4 | Times | Decrease ; 8 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 159000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.55    | mm  | Main Thickness ; | 1.56     | mm           |
| Filling Depth :               | 5.73    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.777   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 10      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 18      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.36 | kN |             |    |
|                  | Max ;     | 0.43 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 41 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 41 |

### Tablet Data (N/A)

|                 |                     |         |       |          |         |
|-----------------|---------------------|---------|-------|----------|---------|
| Total Counter : | 1585878             | Tablets |       |          |         |
| A.W.C Counter : | 1581731             | Tablets |       |          |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 5 (0.0%) | Tablets |
| GOOD :          | 1581726<br>(100.0%) | Tablets |       |          |         |