

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METF 500MG<br>STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2505650             |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 02 Min 18 Sec  |
| Total Counter :  | 1250454       | Total Running Time : | 959 Hour 34 Min        |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 105               | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 180               | Min     | Supply Time ;      | 0.8          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 150               | Min     | Supply Time ;      | 2.2          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 15.66             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 15.05             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 12.28             | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 12.04             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 11.80             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 9.03              | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 8.43              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |                |       |                |       |    |            |    |
|----------------------|----------------|-------|----------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.86           | kN    | Min :          | 2.83  | kN | Punch NO.; | 15 |
| SD :                 | 0.60           | %     | Max :          | 2.91  | kN | Punch NO.; | 11 |
| Mean Main Pressure : | 12.49          | kN    | Min :          | 11.23 | kN | Punch NO.; | 25 |
| SD :                 | 4.70           | %     | Max :          | 13.61 | kN | Punch NO.; | 5  |
| Filling Depth :      | Increase ; 340 | Times | Decrease ; 329 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 29                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.35    | mm  | Pre Thickness ;  | 3.36              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.40    | mm  | Main Thickness ; | 2.41              | mm           |
| Filling Depth :               | 7.64    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.798   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 116     | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 97      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 10 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 24 |
| Ejection Force : | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 1250454             | Tablets |       |           |         |
| A.W.C Counter : | 1242723             | Tablets |       |           |         |
| HEP :           | 249 (0.0%)          | Tablets | LEP : | 36 (0.0%) | Tablets |
| GOOD :          | 1242438<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METF 500MG<br>STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2505650             |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 02 Min 18 Sec  |
| Total Counter :  | 1250454       | Total Running Time : | 959 Hour 34 Min        |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 100               | %       | Manual ;           | 20           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 180               | Min     | Supply Time ;      | 0.8          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 150               | Min     | Supply Time ;      | 2.2          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 15.54             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 14.94             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 12.19             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 11.95             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 11.71             | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 8.96              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 8.37              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.82       | kN  | Min : | 2.78       | kN  | Punch NO.; | 30 |
| SD :                 | 0.70       | %   | Max : | 2.87       | kN  | Punch NO.; | 46 |
| Mean Main Pressure : | 11.87      | kN  | Min : | 10.19      | kN  | Punch NO.; | 39 |
| SD :                 | 6.60       | %   | Max : | 13.74      | kN  | Punch NO.; | 24 |
| Filling Depth :      | Increase ; | 287 | Times | Decrease ; | 293 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 28       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 3.31     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.40    | mm  | Main Thickness ; | 2.41     | mm           |
| Filling Depth :               | 7.42    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.798   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 116     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 97      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 10 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 23 |
| Ejection Force : | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 1250454             | Tablets |       |           |         |
| A.W.C Counter : | 1242723             | Tablets |       |           |         |
| HEP :           | 249 (0.0%)          | Tablets | LEP : | 36 (0.0%) | Tablets |
| GOOD :          | 1242438<br>(100.0%) | Tablets |       |           |         |