

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METF 500MG<br>STEW1A00 |
| User ID :        | ms88737       | Batch NO. :          | EWM2505650             |
| Print Interval : | 50000 Tablets | Running Time :       | 23 Hour 35 Min 11 Sec  |
| Total Counter :  | 4201976       | Total Running Time : | 976 Hour 8 Min         |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 114               | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 120               | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 120               | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 15.27             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 14.68             | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.5 %       | 12.04             | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 11.74             | kN      |                    |              |           |
| LCP :                              | 2.5 %       | 11.45             | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 %      | 8.81              | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 8.22              | kN      | No                 |              |           |
| SD limit :                         | 40.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 35.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |      |       |            |      |            |    |
|----------------------|------------|------|-------|------------|------|------------|----|
| Mean Pre Pressure :  | 3.03       | kN   | Min : | 2.97       | kN   | Punch NO.; | 15 |
| SD :                 | 0.70       | %    | Max : | 3.07       | kN   | Punch NO.; | 40 |
| Mean Main Pressure : | 11.72      | kN   | Min : | 9.98       | kN   | Punch NO.; | 24 |
| SD :                 | 6.50       | %    | Max : | 13.86      | kN   | Punch NO.; | 28 |
| Filling Depth :      | Increase ; | 1170 | Times | Decrease ; | 1118 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 32       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.35    | mm  | Pre Thickness ;  | 3.36     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.40    | mm  | Main Thickness ; | 2.41     | mm           |
| Filling Depth :               | 7.71    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.815   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 78      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 19      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.02 | kN | Max Punch ; | 20 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 37 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 28 |

### Tablet Data (N/A)

|                 |                 |         |       |             |         |
|-----------------|-----------------|---------|-------|-------------|---------|
| Total Counter : | 4201976         | Tablets |       |             |         |
| A.W.C Counter : | 4186175         | Tablets |       |             |         |
| HEP :           | 8293 (0.2%)     | Tablets | LEP : | 1193 (0.0%) | Tablets |
| GOOD :          | 4176689 (99.8%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METF 500MG<br>STEW1A00 |
| User ID :        | ms88737       | Batch NO. :          | EWM2505650             |
| Print Interval : | 50000 Tablets | Running Time :       | 23 Hour 35 Min 11 Sec  |
| Total Counter :  | 4201976       | Total Running Time : | 976 Hour 8 Min         |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 114               | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 4500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 120               | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 120               | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 14.46             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 13.91             | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.5 %  | 11.40             | kN      | Yes 20             | Times        |           |
| Ref :                              |        | 11.12             | kN      |                    |              |           |
| LCP :                              | 2.5 %  | 10.85             | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 % | 8.34              | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 % | 7.79              | kN      | No                 |              |           |
| SD limit :                         | 40.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 35.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |      |       |            |      |            |    |
|----------------------|------------|------|-------|------------|------|------------|----|
| Mean Pre Pressure :  | 2.78       | kN   | Min : | 2.74       | kN   | Punch NO.; | 47 |
| SD :                 | 0.70       | %    | Max : | 2.84       | kN   | Punch NO.; | 22 |
| Mean Main Pressure : | 11.12      | kN   | Min : | 8.52       | kN   | Punch NO.; | 38 |
| SD :                 | 9.60       | %    | Max : | 13.61      | kN   | Punch NO.; | 3  |
| Filling Depth :      | Increase ; | 1314 | Times | Decrease ; | 1325 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 32       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 3.31     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.42    | mm  | Main Thickness ; | 2.43     | mm           |
| Filling Depth :               | 7.19    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.815   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 78      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 19      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 20 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 37 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 28 |

### Tablet Data (N/A)

|                 |                 |         |       |             |         |
|-----------------|-----------------|---------|-------|-------------|---------|
| Total Counter : | 4201976         | Tablets |       |             |         |
| A.W.C Counter : | 4186175         | Tablets |       |             |         |
| HEP :           | 8293 (0.2%)     | Tablets | LEP : | 1193 (0.0%) | Tablets |
| GOOD :          | 4176689 (99.8%) | Tablets |       |             |         |