

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO 400MGSTIU2B00
User ID :	ka98463	Batch NO. :	IUM2500340
Print Interval :	50000 Tablets	Running Time :	03 Hour 23 Min 24 Sec
Total Counter :	500220	Total Running Time :	988 Hour 52 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	200	%	Manual ;	15	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		2025000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	20.66	kN	No		
HEP :	25.0 %	19.87	kN	Yes	15	Rotations
HCP :	2.0 %	16.21	kN	Yes	20	Times
Ref :		15.90	kN			
LCP :	2.0 %	15.58	kN	Yes	20	Times
LEP :	25.0 %	11.92	kN	Yes	15	Rotations
LSP :	30.0 %	11.13	kN	No		
SD limit :	40.0 %			Yes		
Pre HSP		35.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	5.07	kN	Min :	4.61	kN	Punch NO.;	44
		4.00	%	Max :	5.55	kN	Punch NO.;	38
Mean Main Pressure :		15.65	kN	Min :	13.66	kN	Punch NO.;	35
	SD :	5.80	%	Max :	17.51	kN	Punch NO.;	51
Filling Depth :	Increase ;	196	Times	Decrease ;	191	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	159000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	50	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.45	mm	Main Thickness ;	1.46	mm
Filling Depth :	8.29	mm	Current Cam :	Cam C	
Main Air pressure :	0.767	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	13	Min			
Lower Punch Remain Time (S2):	13	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.08	kN		
	Max ;	0.21	kN	Max Punch ;	42
Lower Punch :	Average ;	0.19	kN		
	Max ;	0.64	kN	Max Punch ;	50
Ejection Force :	Average ;	0.06	kN		
	Max ;	0.15	kN	Max Punch ;	15

Tablet Data (N/A)

Total Counter :	500220	Tablets			
A.W.C Counter :	489323	Tablets			
HEP :	74 (0.0%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	489248 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO 400MGSTIU2B00
User ID :	ka98463	Batch NO. :	IUM2500340
Print Interval :	50000 Tablets	Running Time :	03 Hour 23 Min 24 Sec
Total Counter :	500220	Total Running Time :	988 Hour 52 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	200	%	Manual ;	15	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		2025000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	20.54	kN	No		
HEP :	25.0 %	19.75	kN	Yes	15	Rotations
HCP :	2.0 %	16.12	kN	Yes	20	Times
Ref :		15.80	kN			
LCP :	2.0 %	15.49	kN	Yes	20	Times
LEP :	25.0 %	11.85	kN	Yes	15	Rotations
LSP :	30.0 %	11.06	kN	No		
SD limit :	40.0 %			Yes		
Pre HSP		35.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	5.05	kN	Min :	4.64	kN	Punch NO.;	10
		3.70	%	Max :	5.46	kN	Punch NO.;	48
Mean Main Pressure :		15.75	kN	Min :	14.28	kN	Punch NO.;	23
	SD :	5.10	%	Max :	17.64	kN	Punch NO.;	53
Filling Depth :	Increase ;	204	Times	Decrease ;	143	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	159000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	50	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.60	mm	Main Thickness ;	1.61	mm
Filling Depth :	9.01	mm	Current Cam :	Cam C	
Main Air pressure :	0.767	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	13	Min			
Lower Punch Remain Time (S2):	13	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.08	kN		
	Max ;	0.21	kN	Max Punch ;	42
Lower Punch :	Average ;	0.19	kN		
	Max ;	0.65	kN	Max Punch ;	50
Ejection Force :	Average ;	0.06	kN		
	Max ;	0.15	kN	Max Punch ;	15

Tablet Data (N/A)

Total Counter :	500220	Tablets			
A.W.C Counter :	489323	Tablets			
HEP :	74 (0.0%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	489248 (100.0%)	Tablets			