

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO 400MGSTIU2B00
User ID :	am98755	Batch NO. :	IUM2500340
Print Interval :	50000 Tablets	Running Time :	14 Hour 18 Min 41 Sec
Total Counter :	1998230	Total Running Time :	999 Hour 48 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	223	%	Manual ;	15	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		2025000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	22.04	kN	No		
HEP :	25.0 %	21.19	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	17.29	kN	Yes 20	Times	
Ref :		16.95	kN			
LCP :	2.0 %	16.61	kN	Yes 20	Times	
LEP :	25.0 %	12.71	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	11.87	kN	No		
SD limit :	40.0 %			Yes		
Pre HSP		35.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	4.74	kN	Min :	4.47	kN	Punch NO.;	14
SD :	4.70	%	Max :	5.99	kN	Punch NO.;	15
Mean Main Pressure :	13.72	kN	Min :	11.82	kN	Punch NO.;	36
SD :	5.80	%	Max :	15.00	kN	Punch NO.;	19
Filling Depth :	Increase ;	844	Times	Decrease ;	657	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	45	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.45	mm	Main Thickness ;	1.46	mm
Filling Depth :	11.85	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.759	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	12	Min			
Lower Punch Remain Time (S2):	12	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.12	kN		
	Max ;	0.18	kN	Max Punch ;	23
Lower Punch :	Average ;	0.11	kN		
	Max ;	0.69	kN	Max Punch ;	16
Ejection Force :	Average ;	0.07	kN		
	Max ;	0.20	kN	Max Punch ;	16

Tablet Data (N/A)

Total Counter :	1998230	Tablets			
A.W.C Counter :	1979995	Tablets			
HEP :	1469 (0.1%)	Tablets	LEP :	146 (0.0%)	Tablets
GOOD :	1978380 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO 400MGSTIU2B00
User ID :	am98755	Batch NO. :	IUM2500340
Print Interval :	50000 Tablets	Running Time :	14 Hour 18 Min 41 Sec
Total Counter :	1998230	Total Running Time :	999 Hour 48 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	210	%	Manual ;	15	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		2025000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	21.85	kN	No		
HEP :	25.0 %	21.01	kN	Yes	15	Rotations
HCP :	2.0 %	17.14	kN	Yes	20	Times
Ref :		16.80	kN			
LCP :	2.0 %	16.47	kN	Yes	20	Times
LEP :	25.0 %	12.60	kN	Yes	15	Rotations
LSP :	30.0 %	11.76	kN	No		
SD limit :	40.0 %			Yes		
Pre HSP		35.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	4.56	kN	Min :	4.35	kN	Punch NO.;	25
SD :	2.30	%	Max :	4.75	kN	Punch NO.;	50
Mean Main Pressure :	14.47	kN	Min :	13.40	kN	Punch NO.;	32
SD :	4.10	%	Max :	15.88	kN	Punch NO.;	50
Filling Depth :	Increase ;	650	Times	Decrease ;	485	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	42	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.60	mm	Main Thickness ;	1.61	mm
Filling Depth :	10.32	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.759	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	12	Min			
Lower Punch Remain Time (S2):	12	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.13	kN		
	Max ;	0.19	kN	Max Punch ;	23
Lower Punch :	Average ;	0.11	kN		
	Max ;	0.69	kN	Max Punch ;	16
Ejection Force :	Average ;	0.07	kN		
	Max ;	0.20	kN	Max Punch ;	16

Tablet Data (N/A)

Total Counter :	1998230	Tablets			
A.W.C Counter :	1979995	Tablets			
HEP :	1469 (0.1%)	Tablets	LEP :	146 (0.0%)	Tablets
GOOD :	1978380 (99.9%)	Tablets			