

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO 600MGSTIU2C00
User ID :	pb97604	Batch NO. :	IUM2501260
Print Interval :	50000 Tablets	Running Time :	00 Hour 15 Min 27 Sec
Total Counter :	24314	Total Running Time :	1000 Hour 12 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	200	%	Manual ;	19	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		1350000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	23.86	kN	No		
HEP :	25.0 %	22.94	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	18.72	kN	Yes 20	Times	
Ref :		18.36	kN			
LCP :	2.0 %	17.99	kN	Yes 20	Times	
LEP :	25.0 %	13.77	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	12.85	kN	No		
SD limit :	40.0 %			Yes		
Pre HSP		35.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	8.03	kN	Min :	7.25	kN	Punch NO.;	11
		5.50	%	Max :	8.95	kN	Punch NO.;	1
Mean Main Pressure :	SD :	17.92	kN	Min :	16.23	kN	Punch NO.;	15
		5.50	%	Max :	19.88	kN	Punch NO.;	51
Filling Depth :	Increase ;	855	Times	Decrease ;	668	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	30	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	6.30	mm	Pre Thickness ;	4.30	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	5.40	mm	Main Thickness ;	3.40	mm
Filling Depth :	11.56	mm	Current Cam :	Cam C	
Main Air pressure :	0.781	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	15	Min			
Lower Punch Remain Time (S2):	15	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.37	kN		
	Max ;	0.48	kN	Max Punch ;	33
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.11	kN	Max Punch ;	53
Ejection Force :	Average ;	0.05	kN		
	Max ;	0.08	kN	Max Punch ;	15

Tablet Data (N/A)

Total Counter :	24314	Tablets			
A.W.C Counter :	19141	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	19141 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO 600MGSTIU2C00
User ID :	pb97604	Batch NO. :	IUM2501260
Print Interval :	50000 Tablets	Running Time :	00 Hour 15 Min 27 Sec
Total Counter :	24314	Total Running Time :	1000 Hour 12 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	200	%	Manual ;	15	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		1350000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	24.36	kN	No		
HEP :	25.0 %	23.43	kN	Yes	15	Rotations 20 Tablets
HCP :	2.0 %	19.12	kN	Yes	20	Times
Ref :		18.74	kN			
LCP :	2.0 %	18.37	kN	Yes	20	Times
LEP :	25.0 %	14.06	kN	Yes	15	Rotations 20 Tablets
LSP :	30.0 %	13.12	kN	No		
SD limit :	40.0 %			Yes		
Pre HSP		35.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	7.25	kN	Min :	6.65	kN	Punch NO.;	33
SD :	4.20	%	Max :	7.93	kN	Punch NO.;	45
Mean Main Pressure :	18.51	kN	Min :	16.49	kN	Punch NO.;	43
SD :	5.80	%	Max :	20.83	kN	Punch NO.;	51
Filling Depth :	Increase ;	650	Times	Decrease ;	492	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	30	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	6.30	mm	Pre Thickness ;	4.30	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	5.50	mm	Main Thickness ;	3.50	mm
Filling Depth :	10.58	mm	Current Cam :	Cam C	
Main Air pressure :	0.781	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	15	Min			
Lower Punch Remain Time (S2):	15	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.37	kN		
	Max ;	0.48	kN	Max Punch ;	33
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.10	kN	Max Punch ;	53
Ejection Force :	Average ;	0.05	kN		
	Max ;	0.08	kN	Max Punch ;	15

Tablet Data (N/A)

Total Counter :	24314	Tablets			
A.W.C Counter :	19141	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	19141 (100.0%)	Tablets			