

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO 600MGSTIU2C00
User ID :	pb97604	Batch NO. :	IUM2501260
Print Interval :	50000 Tablets	Running Time :	00 Hour 31 Min 46 Sec
Total Counter :	50014	Total Running Time :	1000 Hour 29 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	255	%	Manual ;	19	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		1350000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	24.78	kN	No		
HEP :	25.0 %	23.82	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	19.44	kN	Yes 20	Times	
Ref :		19.06	kN			
LCP :	2.0 %	18.68	kN	Yes 20	Times	
LEP :	25.0 %	14.29	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	13.34	kN	No		
SD limit :	40.0 %			Yes		
Pre HSP		35.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	7.91	kN	Min :	6.97	kN	Punch NO.;	30
SD :	5.40	%	Max :	8.66	kN	Punch NO.;	15
Mean Main Pressure :	18.69	kN	Min :	16.49	kN	Punch NO.;	44
SD :	5.10	%	Max :	20.97	kN	Punch NO.;	51
Filling Depth :	Increase ;	872	Times	Decrease ;	686	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	38	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	6.30	mm	Pre Thickness ;	4.30	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	5.35	mm	Main Thickness ;	3.35	mm
Filling Depth :	11.56	mm	Current Cam :	Cam C	
Main Air pressure :	0.781	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	14	Min			
Lower Punch Remain Time (S2):	14	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.36	kN		
	Max ;	0.48	kN	Max Punch ;	27
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.13	kN	Max Punch ;	15
Ejection Force :	Average ;	0.05	kN		
	Max ;	0.08	kN	Max Punch ;	15

Tablet Data (N/A)

Total Counter :	50014	Tablets			
A.W.C Counter :	44846	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	44846 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO 600MGSTIU2C00
User ID :	pb97604	Batch NO. :	IUM2501260
Print Interval :	50000 Tablets	Running Time :	00 Hour 31 Min 46 Sec
Total Counter :	50014	Total Running Time :	1000 Hour 29 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	211	%	Manual ;	15	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		1350000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	26.00	kN	No		
HEP :	25.0 %	25.00	kN	Yes	15	Rotations
HCP :	2.0 %	20.40	kN	Yes	20	Times
Ref :		20.00	kN			
LCP :	2.0 %	19.60	kN	Yes	20	Times
LEP :	25.0 %	15.00	kN	Yes	15	Rotations
LSP :	30.0 %	14.00	kN	No		
SD limit :	40.0 %			Yes		
Pre HSP		35.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	7.81	kN	Min :	7.23	kN	Punch NO.;	20
SD :	3.80	%	Max :	8.58	kN	Punch NO.;	4
Mean Main Pressure :	20.03	kN	Min :	18.19	kN	Punch NO.;	36
SD :	5.00	%	Max :	22.37	kN	Punch NO.;	51
Filling Depth :	Increase ;	654	Times	Decrease ;	494	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	32	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	6.30	mm	Pre Thickness ;	4.30	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	5.45	mm	Main Thickness ;	3.45	mm
Filling Depth :	10.74	mm	Current Cam :	Cam C	
Main Air pressure :	0.781	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	14	Min			
Lower Punch Remain Time (S2):	14	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.37	kN		
	Max ;	0.48	kN	Max Punch ;	27
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.13	kN	Max Punch ;	15
Ejection Force :	Average ;	0.05	kN		
	Max ;	0.08	kN	Max Punch ;	15

Tablet Data (N/A)

Total Counter :	50014	Tablets			
A.W.C Counter :	44846	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	44846 (100.0%)	Tablets			