

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBUPRO<br>600MGSTIU2C00 |
| User ID :        | rp98927       | Batch NO. :          | IUM2501360              |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 48 Min 08 Sec   |
| Total Counter :  | 650218        | Total Running Time : | 1017 Hour 52 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 495               | %       | Manual ;           | 19           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 1350000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 37.48             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 36.04             | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 29.41             | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 28.83             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 28.26             | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 %      | 21.62             | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 20.18             | kN      | No                 |              |           |
| SD limit :                         | 40.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 35.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 8.23       | kN  | Min : | 7.65       | kN  | Punch NO.; | 33 |
| SD :                 | 3.90       | %   | Max : | 8.86       | kN  | Punch NO.; | 42 |
| Mean Main Pressure : | 29.10      | kN  | Min : | 26.80      | kN  | Punch NO.; | 31 |
| SD :                 | 4.50       | %   | Max : | 31.89      | kN  | Punch NO.; | 42 |
| Filling Depth :      | Increase ; | 547 | Times | Decrease ; | 460 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 18      | RPM | Capacity :       | 114480   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 89       | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 4.30     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.35    | mm  | Main Thickness ; | 3.35     | mm           |
| Filling Depth :               | 11.00   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.799   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 11      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 14      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.20 | kN | Max Punch ; | 15 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 16 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 650218            | Tablets |       |            |         |
| A.W.C Counter : | 643881            | Tablets |       |            |         |
| HEP :           | 53 (0.0%)         | Tablets | LEP : | 432 (0.1%) | Tablets |
| GOOD :          | 643396<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBUPRO<br>600MGSTIU2C00 |
| User ID :        | rp98927       | Batch NO. :          | IUM2501360              |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 48 Min 08 Sec   |
| Total Counter :  | 650218        | Total Running Time : | 1017 Hour 52 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 445               | %       | Manual ;           | 50           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 1350000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 60                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 38.75             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 37.26             | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %  | 30.41             | kN      | Yes 20             | Times        |           |
| Ref :                              |        | 29.81             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 29.21             | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 % | 22.36             | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 % | 20.87             | kN      | No                 |              |           |
| SD limit :                         | 40.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 35.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 8.04       | kN  | Min : | 7.47       | kN  | Punch NO.; | 21 |
| SD :                 | 3.30       | %   | Max : | 8.58       | kN  | Punch NO.; | 42 |
| Mean Main Pressure : | 30.13      | kN  | Min : | 27.85      | kN  | Punch NO.; | 14 |
| SD :                 | 4.20       | %   | Max : | 32.80      | kN  | Punch NO.; | 4  |
| Filling Depth :      | Increase ; | 277 | Times | Decrease ; | 220 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 18      | RPM | Capacity :       | 114480   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 80       | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 4.30     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.45    | mm  | Main Thickness ; | 3.45     | mm           |
| Filling Depth :               | 10.58   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.799   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 11      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 14      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.20 | kN | Max Punch ; | 15 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 16 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 650218            | Tablets |       |            |         |
| A.W.C Counter : | 643881            | Tablets |       |            |         |
| HEP :           | 53 (0.0%)         | Tablets | LEP : | 432 (0.1%) | Tablets |
| GOOD :          | 643396<br>(99.9%) | Tablets |       |            |         |