

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Sertra 100mgSTSH4B00  |
| User ID :        | pn98247       | Batch NO. :          | SHM2501310            |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 52 Min 47 Sec |
| Total Counter :  | 2301616       | Total Running Time : | 1044 Hour 43 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 26       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 3000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.500    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 45       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 20.98    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 20.18    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 16.46    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 16.14    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 15.82    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 12.11    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 11.30    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 4.48           | kN    | Min :      | 4.34  | kN    | Punch NO.; | 25 |
| SD :                 | 1.20           | %     | Max :      | 4.60  | kN    | Punch NO.; | 76 |
| Mean Main Pressure : | 16.10          | kN    | Min :      | 15.10 | kN    | Punch NO.; | 25 |
| SD :                 | 2.90           | %     | Max :      | 17.49 | kN    | Punch NO.; | 76 |
| Filling Depth :      | Increase ; 249 | Times | Decrease ; | 302   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 7        | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.60    | mm  | Pre Thickness ;  | 3.10     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.53    | mm  | Main Thickness ; | 2.03     | mm           |
| Filling Depth :               | 7.80    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.727   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 32      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.23 | kN |             |    |
|                  | Max ;     | 0.36 | kN | Max Punch ; | 7  |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 37 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 74 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 2301616             | Tablets |       |           |         |
| A.W.C Counter : | 2298254             | Tablets |       |           |         |
| HEP :           | 8 (0.0%)            | Tablets | LEP : | 22 (0.0%) | Tablets |
| GOOD :          | 2298224<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Sertra 100mgSTSH4B00  |
| User ID :        | pn98247       | Batch NO. :          | SHM2501310            |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 52 Min 47 Sec |
| Total Counter :  | 2301616       | Total Running Time : | 1044 Hour 43 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 25       | %       | Manual ;           | 48           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 3000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.500    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 45       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 21.47    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 20.64    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 16.84    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 16.51    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 16.18    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 12.38    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 11.56    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |      |       |            |      |            |    |
|----------------------|------------|------|-------|------------|------|------------|----|
| Mean Pre Pressure :  | 4.53       | kN   | Min : | 4.44       | kN   | Punch NO.; | 19 |
| SD :                 | 1.20       | %    | Max : | 4.71       | kN   | Punch NO.; | 76 |
| Mean Main Pressure : | 15.87      | kN   | Min : | 14.94      | kN   | Punch NO.; | 12 |
| SD :                 | 3.90       | %    | Max : | 18.27      | kN   | Punch NO.; | 76 |
| Filling Depth :      | Increase ; | 1065 | Times | Decrease ; | 1056 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 7        | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.60    | mm  | Pre Thickness ;  | 3.10     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.64    | mm  | Main Thickness ; | 2.14     | mm           |
| Filling Depth :               | 7.74    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.727   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 32      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.22 | kN |             |    |
|                  | Max ;     | 0.36 | kN | Max Punch ; | 7  |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 37 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 74 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 2301616             | Tablets |       |           |         |
| A.W.C Counter : | 2298254             | Tablets |       |           |         |
| HEP :           | 8 (0.0%)            | Tablets | LEP : | 22 (0.0%) | Tablets |
| GOOD :          | 2298224<br>(100.0%) | Tablets |       |           |         |