

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MET 500 MG,STEW1A00
User ID :	am98755	Batch NO. :	EWM2507950
Print Interval :	50000 Tablets	Running Time :	00 Hour 21 Min 46 Sec
Total Counter :	52102	Total Running Time :	1109 Hour 35 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	60	%	Manual ;	15	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		4500000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		45	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		18	Min	Supply Time ;	3.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	13.26	kN	No		
HEP :	25.0 %	12.75	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	10.40	kN	Yes	15	Times
Ref :		10.20	kN			
LCP :	2.0 %	9.99	kN	Yes	15	Times
LEP :	25.0 %	7.65	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	7.14	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.78	kN	Min :	3.57	kN	Punch NO.;	29
SD :	3.00	%	Max :	4.07	kN	Punch NO.;	9
Mean Main Pressure :	10.38	kN	Min :	9.28	kN	Punch NO.;	3
SD :	5.00	%	Max :	11.62	kN	Punch NO.;	21
Filling Depth :	Increase ;	17	Times	Decrease ;	5	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	159000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	15	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	6.30	mm	Pre Thickness ;	3.31	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	5.05	mm	Main Thickness ;	2.06	mm
Filling Depth :	7.21	mm	Current Cam :	Cam B	
Main Air pressure :	0.808	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	7	Min			
Lower Punch Remain Time (S2):	17	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.06	kN		
	Max ;	0.22	kN	Max Punch ;	41
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.20	kN	Max Punch ;	7
Ejection Force :	Average ;	0.06	kN		
	Max ;	0.15	kN	Max Punch ;	53

Tablet Data (N/A)

Total Counter :	52102	Tablets			
A.W.C Counter :	40695	Tablets			
HEP :	5 (0.0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	40690 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MET 500 MG,STEW1A00
User ID :	am98755	Batch NO. :	EWM2507950
Print Interval :	50000 Tablets	Running Time :	00 Hour 21 Min 46 Sec
Total Counter :	52102	Total Running Time :	1109 Hour 35 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	60	%	Manual ;	20	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		4500000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		45	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		18	Min	Supply Time ;	3.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	12.99	kN	No		
HEP :	25.0 %	12.49	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	10.19	kN	Yes	15	Times
Ref :		9.99	kN			
LCP :	2.0 %	9.79	kN	Yes	15	Times
LEP :	25.0 %	7.49	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	6.99	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	3.38	kN	Min :	3.09	kN	Punch NO.;	26
		3.50	%	Max :	3.65	kN	Punch NO.;	51
Mean Main Pressure :	SD :	10.21	kN	Min :	9.23	kN	Punch NO.;	26
		5.20	%	Max :	11.71	kN	Punch NO.;	21
Filling Depth :	Increase ;	6	Times	Decrease ;	6	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	159000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	15	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	6.30	mm	Pre Thickness ;	3.31	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	5.25	mm	Main Thickness ;	2.26	mm
Filling Depth :	6.58	mm	Current Cam :	Cam B	
Main Air pressure :	0.808	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	7	Min			
Lower Punch Remain Time (S2):	17	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.06	kN		
	Max ;	0.21	kN	Max Punch ;	45
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.20	kN	Max Punch ;	7
Ejection Force :	Average ;	0.06	kN		
	Max ;	0.15	kN	Max Punch ;	53

Tablet Data (N/A)

Total Counter :	52102	Tablets			
A.W.C Counter :	40695	Tablets			
HEP :	5 (0.0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	40690 (100.0%)	Tablets			