

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2508350             |
| Print Interval : | 50000 Tablets | Running Time :       | 22 Hour 57 Min 33 Sec  |
| Total Counter :  | 4095738       | Total Running Time : | 1190 Hour 30 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 91                | %       | Manual ;           | 15           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 60                | Min     | Supply Time ;      | 2.0          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 100               | Min     | Supply Time ;      | 3.0          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 10.34             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 9.94              | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 8.11              | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 7.96              | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 7.80              | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 5.97              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 5.57              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | None        |                   |         |                    |              |            |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.96       | kN | Min : | 2.81       | kN | Punch NO.; | 16 |
| SD :                 | 2.50       | %  | Max : | 3.17       | kN | Punch NO.; | 8  |
| Mean Main Pressure : | 8.04       | kN | Min : | 7.21       | kN | Punch NO.; | 43 |
| SD :                 | 4.70       | %  | Max : | 8.97       | kN | Punch NO.; | 21 |
| Filling Depth :      | Increase ; | 88 | Times | Decrease ; | 98 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 25       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.46    | mm  | Pre Thickness ;  | 3.46     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.26    | mm  | Main Thickness ; | 2.26     | mm           |
| Filling Depth :               | 7.30    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.849   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 33      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 94      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 29 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 14 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 4095738            | Tablets |       |            |         |
| A.W.C Counter : | 4087893            | Tablets |       |            |         |
| HEP :           | 2102 (0.1%)        | Tablets | LEP : | 722 (0.0%) | Tablets |
| GOOD :          | 4085069<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2508350             |
| Print Interval : | 50000 Tablets | Running Time :       | 22 Hour 57 Min 33 Sec  |
| Total Counter :  | 4095738       | Total Running Time : | 1190 Hour 30 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 91                | %       | Manual ;           | 20           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 60                | Min     | Supply Time ;      | 2.0          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 100               | Min     | Supply Time ;      | 3.0          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 10.33             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 9.93              | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 8.10              | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 7.95              | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 7.79              | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 5.96              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 5.56              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.90       | kN | Min : | 2.73       | kN | Punch NO.; | 13 |
| SD :                 | 2.70       | %  | Max : | 3.11       | kN | Punch NO.; | 47 |
| Mean Main Pressure : | 7.94       | kN | Min : | 7.08       | kN | Punch NO.; | 43 |
| SD :                 | 5.80       | %  | Max : | 8.94       | kN | Punch NO.; | 39 |
| Filling Depth :      | Increase ; | 72 | Times | Decrease ; | 74 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 25       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.39    | mm  | Pre Thickness ;  | 3.40     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.46    | mm  | Main Thickness ; | 2.47     | mm           |
| Filling Depth :               | 6.91    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.849   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 33      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 94      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 29 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.20 | kN | Max Punch ; | 14 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 4095738            | Tablets |       |            |         |
| A.W.C Counter : | 4087893            | Tablets |       |            |         |
| HEP :           | 2102 (0.1%)        | Tablets | LEP : | 722 (0.0%) | Tablets |
| GOOD :          | 4085069<br>(99.9%) | Tablets |       |            |         |