

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | vg96371       | Batch NO. :          | EWM2508950             |
| Print Interval : | 50000 Tablets | Running Time :       | 14 Hour 47 Min 17 Sec  |
| Total Counter :  | 2552086       | Total Running Time : | 1207 Hour 30 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 85                | %       | Manual ;           | 15           | RPM        |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 12.00 mm   |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 120               | Min     | Supply Time ;      | 1.5          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 120               | Min     | Supply Time ;      | 1.8          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 30.45             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 29.28             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 23.89             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 23.42             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 22.95             | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 17.57             | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 16.40             | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | None        |                   |         |                    |              |            |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.89  | kN    | Min :      | 2.70  | kN    | Punch NO.; | 17 |
|                      |            | 2.70  | %     | Max :      | 3.08  | kN    | Punch NO.; | 51 |
| Mean Main Pressure : | SD :       | 23.18 | kN    | Min :      | 21.43 | kN    | Punch NO.; | 29 |
|                      |            | 3.40  | %     | Max :      | 24.88 | kN    | Punch NO.; | 26 |
| Filling Depth :      | Increase ; | 77    | Times | Decrease ; | 67    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 19       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.46    | mm  | Pre Thickness ;  | 3.46     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.20    | mm  | Main Thickness ; | 2.20     | mm           |
| Filling Depth :               | 7.45    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.838   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 58      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 59      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 2552086             | Tablets |       |            |         |
| A.W.C Counter : | 2536561             | Tablets |       |            |         |
| HEP :           | 154 (0.0%)          | Tablets | LEP : | 310 (0.0%) | Tablets |
| GOOD :          | 2536097<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | vg96371       | Batch NO. :          | EWM2508950             |
| Print Interval : | 50000 Tablets | Running Time :       | 14 Hour 47 Min 17 Sec  |
| Total Counter :  | 2552086       | Total Running Time : | 1207 Hour 30 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 80        | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |        | 8.5-14mm  |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |        | 4500000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.000     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 120       | Min     | Supply Time ;      | 1.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 120       | Min     | Supply Time ;      | 1.8          | Sec       |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 27.16     | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 26.12     | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 21.31     | kN      | Yes 10             | Times        |           |
| Ref :                              |        | 20.89     | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 20.48     | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 % | 15.67     | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 14.63     | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 30.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.64       | kN | Min : | 2.40       | kN | Punch NO.; | 27 |
| SD :                 | 4.20       | %  | Max : | 2.86       | kN | Punch NO.; | 40 |
| Mean Main Pressure : | 20.57      | kN | Min : | 18.57      | kN | Punch NO.; | 29 |
| SD :                 | 3.90       | %  | Max : | 21.97      | kN | Punch NO.; | 1  |
| Filling Depth :      | Increase ; | 56 | Times | Decrease ; | 59 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.40    | mm  | Pre Thickness ;  | 3.41     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.40    | mm  | Main Thickness ; | 2.41     | mm           |
| Filling Depth :               | 7.08    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.838   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 58      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 59      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 2552086             | Tablets |       |            |         |
| A.W.C Counter : | 2536561             | Tablets |       |            |         |
| HEP :           | 154 (0.0%)          | Tablets | LEP : | 310 (0.0%) | Tablets |
| GOOD :          | 2536097<br>(100.0%) | Tablets |       |            |         |