

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | sj97904       | Batch NO. :          | EWM2509450             |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 13 Min 56 Sec  |
| Total Counter :  | 1467368       | Total Running Time : | 1227 Hour 9 Min        |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 68                | %       | Manual ;           | 15           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 20                | Min     | Supply Time ;      | 1.0          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 1.0          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 13.33             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 12.81             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 10.46             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 10.25             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 10.05             | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 7.69              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 7.18              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | None        |                   |         |                    |              |            |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.91       | kN  | Min : | 2.77       | kN  | Punch NO.; | 16 |
| SD :                 | 2.40       | %   | Max : | 3.08       | kN  | Punch NO.; | 30 |
| Mean Main Pressure : | 10.02      | kN  | Min : | 8.59       | kN  | Punch NO.; | 30 |
| SD :                 | 5.80       | %   | Max : | 11.18      | kN  | Punch NO.; | 26 |
| Filling Depth :      | Increase ; | 413 | Times | Decrease ; | 405 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 19       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.46    | mm  | Pre Thickness ;  | 3.46     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.34    | mm  | Main Thickness ; | 2.35     | mm           |
| Filling Depth :               | 7.08    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.806   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 1       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 1       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.28 | kN | Max Punch ; | 42 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.22 | kN |             |    |
|                  | Max ;     | 0.39 | kN | Max Punch ; | 10 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1467368            | Tablets |       |            |         |
| A.W.C Counter : | 1461609            | Tablets |       |            |         |
| HEP :           | 510 (0.0%)         | Tablets | LEP : | 442 (0.0%) | Tablets |
| GOOD :          | 1460657<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | sj97904       | Batch NO. :          | EWM2509450             |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 13 Min 56 Sec  |
| Total Counter :  | 1467368       | Total Running Time : | 1227 Hour 9 Min        |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 64                | %       | Manual ;           | 20           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 20                | Min     | Supply Time ;      | 1.0          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 1.0          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 12.74             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 12.25             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 10.00             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 9.80              | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 9.61              | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 7.35              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 6.86              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2 Rotations  |            |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.65       | kN  | Min : | 2.49       | kN  | Punch NO.; | 27 |
| SD :                 | 2.60       | %   | Max : | 2.81       | kN  | Punch NO.; | 40 |
| Mean Main Pressure : | 9.47       | kN  | Min : | 7.97       | kN  | Punch NO.; | 29 |
| SD :                 | 6.30       | %   | Max : | 10.79      | kN  | Punch NO.; | 20 |
| Filling Depth :      | Increase ; | 166 | Times | Decrease ; | 152 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.40    | mm  | Pre Thickness ;  | 3.41     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.48    | mm  | Main Thickness ; | 2.49     | mm           |
| Filling Depth :               | 6.91    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.806   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 1       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 1       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.28 | kN | Max Punch ; | 42 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.22 | kN |             |    |
|                  | Max ;     | 0.36 | kN | Max Punch ; | 24 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1467368            | Tablets |       |            |         |
| A.W.C Counter : | 1461609            | Tablets |       |            |         |
| HEP :           | 510 (0.0%)         | Tablets | LEP : | 442 (0.0%) | Tablets |
| GOOD :          | 1460657<br>(99.9%) | Tablets |       |            |         |