

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA50MG,STJL1A0<br>0 |
| User ID :        | ks99394       | Batch NO. :          | JLM2500450              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 24 Min 21 Sec   |
| Total Counter :  | 55876         | Total Running Time : | 1551 Hour 44 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 35                | %       | Manual ;           | 35           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 700000            | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.500             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 70                | Min     | Supply Time ;      | 0.2          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 60                | Min     | Supply Time ;      | 0.2          | Sec        |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 8.33              | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 8.01              | kN      | Yes                | 10 Rotations | 15 Tablets |
| HCP :                              | 2.0 %       | 6.54              | kN      | Yes                | 15 Times     |            |
| Ref :                              |             | 6.41              | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 6.28              | kN      | Yes                | 15 Times     |            |
| LEP :                              | 25.0 %      | 4.81              | kN      | Yes                | 10 Rotations | 15 Tablets |
| LSP :                              | 30.0 %      | 4.49              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations  |
| Empty Punch No. :                  | None        |                   |         |                    |              |            |

### Pressure Data

|                      |              |       |              |       |    |            |    |
|----------------------|--------------|-------|--------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.60         | kN    | Min :        | 2.58  | kN | Punch NO.; | 17 |
| SD :                 | 0.30         | %     | Max :        | 2.63  | kN | Punch NO.; | 76 |
| Mean Main Pressure : | 6.47         | kN    | Min :        | 6.04  | kN | Punch NO.; | 16 |
| SD :                 | 2.70         | %     | Max :        | 7.14  | kN | Punch NO.; | 76 |
| Filling Depth :      | Increase ; 4 | Times | Decrease ; 7 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 138600   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 5        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.49    | mm  | Pre Thickness ;  | 3.50     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.70    | mm  | Main Thickness ; | 1.71     | mm           |
| Filling Depth :               | 5.47    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.761   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 44      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.42 | kN |             |    |
|                  | Max ;     | 0.68 | kN | Max Punch ; | 12 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 61 |
| Ejection Force : | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 1  |

### Tablet Data (N/A)

|                 |               |         |       |           |         |
|-----------------|---------------|---------|-------|-----------|---------|
| Total Counter : | 55876         | Tablets |       |           |         |
| A.W.C Counter : | 49129         | Tablets |       |           |         |
| HEP :           | 215 (0.4%)    | Tablets | LEP : | 60 (0.1%) | Tablets |
| GOOD :          | 48854 (99.4%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA50MG,STJL1A00   |
| User ID :        | ks99394       | Batch NO. :          | JLM2500450            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 24 Min 21 Sec |
| Total Counter :  | 55876         | Total Running Time : | 1551 Hour 44 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 35        | %       | Manual ;           | 16           | RPM       |
| Filling Cam :                      |        | Cam B     |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |        | 4.5-10mm  |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |        | 700000    | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.500     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 70        | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 60        | Min     | Supply Time ;      | 0.2          | Sec       |
| Powder Supply Time :               |        | 30.0      | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 5.91      | kN      | No                 |              |           |
| HEP :                              | 24.0 % | 5.63      | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 4.63      | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 4.54      | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 4.45      | kN      | Yes 15             | Times        |           |
| LEP :                              | 24.0 % | 3.45      | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 3.18      | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 30.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.37       | kN | Min : | 2.35       | kN | Punch NO.; | 10 |
| SD :                 | 0.50       | %  | Max : | 2.42       | kN | Punch NO.; | 36 |
| Mean Main Pressure : | 4.46       | kN | Min : | 3.97       | kN | Punch NO.; | 19 |
| SD :                 | 5.90       | %  | Max : | 5.61       | kN | Punch NO.; | 76 |
| Filling Depth :      | Increase ; | 11 | Times | Decrease ; | 11 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 138600   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 5        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.49    | mm  | Pre Thickness ;  | 3.50     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.79    | mm  | Main Thickness ; | 1.80     | mm           |
| Filling Depth :               | 5.80    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.761   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 44      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.42 | kN |             |    |
|                  | Max ;     | 0.68 | kN | Max Punch ; | 12 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 61 |
| Ejection Force : | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 25 |

### Tablet Data (N/A)

|                 |                  |         |       |           |         |
|-----------------|------------------|---------|-------|-----------|---------|
| Total Counter : | 55876            | Tablets |       |           |         |
| A.W.C Counter : | 49129            | Tablets |       |           |         |
| HEP :           | 215 (0.4%)       | Tablets | LEP : | 60 (0.1%) | Tablets |
| GOOD :          | 48854<br>(99.4%) | Tablets |       |           |         |