

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBU 400 MG STIR1B00   |
| User ID :        | pb97604       | Batch NO. :          | IRM2507840            |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 09 Min 15 Sec |
| Total Counter :  | 356212        | Total Running Time : | 1571 Hour 46 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |       |           |
|------------------------------------|--------|-----------|---------|--------------------|-------|-----------|
| Feeder :                           | Auto ; | 144       | %       | Manual ;           | 33    | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min : | 6.50 mm   |
| Target Quantity :                  |        | 8.5-14mm  |         |                    | Max : | 12.00 mm  |
| Air Pressure Low Limit :           |        | 1000000   | Tablets |                    |       |           |
| Oil Lubrication :                  |        | 0.400     | MPa     |                    |       |           |
| Upper Punch Interval Time(S1) ;    |        | 25        | Min     | Supply Time ;      | 1.0   | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 15        | Min     | Supply Time ;      | 2.0   | Sec       |
| Powder Supply Time :               |        | 30.0      | Sec     |                    |       |           |
| Initial Reject Time :              |        | 1         | Tablets |                    |       |           |
| Upper Punch Tightness High Limit : |        | 2.00      | kN      |                    |       |           |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |       |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |       |           |
|                                    |        |           |         | Stop Condition ;   |       |           |
| HSP :                              | 30.0 % | 18.76     | kN      | No                 |       |           |
| HEP :                              | 20.0 % | 17.32     | kN      | Yes                | 10    | Rotations |
| HCP :                              | 2.0 %  | 14.72     | kN      | Yes                | 15    | Times     |
| Ref :                              |        | 14.43     | kN      |                    |       |           |
| LCP :                              | 2.0 %  | 14.14     | kN      | Yes                | 15    | Times     |
| LEP :                              | 20.0 % | 11.55     | kN      | Yes                | 10    | Rotations |
| LSP :                              | 30.0 % | 10.10     | kN      | No                 |       |           |
| SD limit :                         | 35.0 % |           |         | Yes                |       |           |
| Pre HSP                            |        | 30.00     | kN      | Yes                |       |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 1     | Rotations |
| Empty Punch No. :                  | None   |           |         |                    |       |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.63  | kN    | Min :      | 1.48  | kN    | Punch NO.; | 2  |
|                      |            | 4.40  | %     | Max :      | 1.74  | kN    | Punch NO.; | 48 |
| Mean Main Pressure : |            | 14.43 | kN    | Min :      | 14.36 | kN    | Punch NO.; | 22 |
|                      | SD :       | 0.30  | %     | Max :      | 14.58 | kN    | Punch NO.; | 26 |
| Filling Depth :      | Increase ; | 0     | Times | Decrease ; | 0     | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 18     | RPM | Capacity :       | 114480   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 26       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.20   | mm  | Pre Thickness ;  | 2.21     | mm           |
| Main Upper Punch Position ;   | 3.00   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.42   | mm  | Main Thickness ; | 1.42     | mm           |
| Filling Depth :               | 7.49   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.725  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 9      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 10     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.53 | kN |             |    |
|                  | Max ;     | 0.60 | kN | Max Punch ; | 11 |
| Lower Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.03 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.34 | kN | Max Punch ; | 17 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 356212             | Tablets |       |        |         |
| A.W.C Counter : | 347887             | Tablets |       |        |         |
| HEP :           | 130 (0.0%)         | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 347757<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBU 400 MG STIR1B00   |
| User ID :        | pb97604       | Batch NO. :          | IRM2507840            |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 09 Min 15 Sec |
| Total Counter :  | 356212        | Total Running Time : | 1571 Hour 46 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 145      | %       | Manual ;           | 16           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 1000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 30.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 18.68    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 17.24    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 14.66    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 14.37    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 14.08    | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 %      | 11.49    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 10.06    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |              |       |              |       |    |            |    |
|----------------------|--------------|-------|--------------|-------|----|------------|----|
| Mean Pre Pressure :  | 1.89         | kN    | Min :        | 1.72  | kN | Punch NO.; | 17 |
| SD :                 | 4.50         | %     | Max :        | 2.17  | kN | Punch NO.; | 36 |
| Mean Main Pressure : | 14.37        | kN    | Min :        | 14.29 | kN | Punch NO.; | 53 |
| SD :                 | 0.20         | %     | Max :        | 14.45 | kN | Punch NO.; | 36 |
| Filling Depth :      | Increase ; 0 | Times | Decrease ; 0 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 18     | RPM | Capacity :       | 114480   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 26       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.20   | mm  | Pre Thickness ;  | 2.21     | mm           |
| Main Upper Punch Position ;   | 3.00   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.45   | mm  | Main Thickness ; | 1.45     | mm           |
| Filling Depth :               | 7.69   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.725  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 9      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 10     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.52 | kN |             |    |
|                  | Max ;     | 0.61 | kN | Max Punch ; | 11 |
| Lower Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.03 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.34 | kN | Max Punch ; | 17 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 356212             | Tablets |       |        |         |
| A.W.C Counter : | 347887             | Tablets |       |        |         |
| HEP :           | 130 (0.0%)         | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 347757<br>(100.0%) | Tablets |       |        |         |