

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBU 400 MG STIR1B00   |
| User ID :        | pb97604       | Batch NO. :          | IRM2507840            |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 53 Min 52 Sec |
| Total Counter :  | 656314        | Total Running Time : | 1574 Hour 31 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 186      | %       | Manual ;           | 33           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 1000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 30.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 18.75    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 17.31    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 14.71    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 14.43    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 14.14    | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 %      | 11.54    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 10.10    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.66  | kN    | Min :      | 1.50  | kN    | Punch NO.; | 45 |
|                      |            | 4.90  | %     | Max :      | 1.87  | kN    | Punch NO.; | 11 |
| Mean Main Pressure : |            | 14.42 | kN    | Min :      | 14.36 | kN    | Punch NO.; | 2  |
|                      | SD :       | 0.20  | %     | Max :      | 14.52 | kN    | Punch NO.; | 24 |
| Filling Depth :      | Increase ; | 0     | Times | Decrease ; | 0     | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 16     | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 30                | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.20   | mm  | Pre Thickness ;  | 2.21              | mm           |
| Main Upper Punch Position ;   | 3.00   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.42   | mm  | Main Thickness ; | 1.42              | mm           |
| Filling Depth :               | 7.39   | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.759  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 22     | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 13     | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.36 | kN |             |    |
|                  | Max ;     | 0.40 | kN | Max Punch ; | 37 |
| Lower Punch :    | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.03 | kN | Max Punch ; | 1  |
| Ejection Force : | Average ; | 0.16 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 33 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 656314             | Tablets |       |        |         |
| A.W.C Counter : | 647921             | Tablets |       |        |         |
| HEP :           | 130 (0.0%)         | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 647791<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBU 400 MG STIR1B00   |
| User ID :        | pb97604       | Batch NO. :          | IRM2507840            |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 53 Min 52 Sec |
| Total Counter :  | 656314        | Total Running Time : | 1574 Hour 31 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 176      | %       | Manual ;           | 16           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 1000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 30.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 18.71    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 17.27    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 14.68    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 14.39    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 14.10    | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 %      | 11.51    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 10.07    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.63  | kN    | Min :      | 1.43  | kN    | Punch NO.; | 17 |
|                      |            | 5.70  | %     | Max :      | 1.88  | kN    | Punch NO.; | 50 |
| Mean Main Pressure : |            | 14.38 | kN    | Min :      | 14.33 | kN    | Punch NO.; | 51 |
|                      | SD :       | 0.20  | %     | Max :      | 14.46 | kN    | Punch NO.; | 24 |
| Filling Depth :      | Increase ; | 0     | Times | Decrease ; | 0     | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 16     | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 28                | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.20   | mm  | Pre Thickness ;  | 2.21              | mm           |
| Main Upper Punch Position ;   | 3.00   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.45   | mm  | Main Thickness ; | 1.45              | mm           |
| Filling Depth :               | 7.69   | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.759  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 22     | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 13     | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.36 | kN |             |    |
|                  | Max ;     | 0.40 | kN | Max Punch ; | 37 |
| Lower Punch :    | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.03 | kN | Max Punch ; | 1  |
| Ejection Force : | Average ; | 0.16 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 33 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 656314             | Tablets |       |        |         |
| A.W.C Counter : | 647921             | Tablets |       |        |         |
| HEP :           | 130 (0.0%)         | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 647791<br>(100.0%) | Tablets |       |        |         |