

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 20 MG,STKQ1A00
User ID :	rp98927	Batch NO. :	KQM2501920
Print Interval :	50000 Tablets	Running Time :	00 Hour 39 Min 27 Sec
Total Counter :	150086	Total Running Time :	1578 Hour 31 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	15	%	Manual ;	18	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		5500000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		25	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		30.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	12.20	kN	No		
HEP :	20.0 %	11.26	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	9.57	kN	Yes 15	Times	
Ref :		9.38	kN			
LCP :	2.0 %	9.19	kN	Yes 15	Times	
LEP :	20.0 %	7.51	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	6.57	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	1.29	kN	Min :	1.03	kN	Punch NO.;	16
SD :	10.60	%	Max :	1.68	kN	Punch NO.;	54
Mean Main Pressure :	8.65	kN	Min :	7.63	kN	Punch NO.;	65
SD :	5.20	%	Max :	9.78	kN	Punch NO.;	54
Filling Depth :	Increase ;	151	Times	Decrease ;	150	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	30	RPM	Capacity :	234000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	4	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.10	mm	Pre Thickness ;	1.10	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.55	mm	Main Thickness ;	0.55	mm
Filling Depth :	6.66	mm	Current Cam :	Cam B	
Main Air pressure :	0.718	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	21	Min			
Lower Punch Remain Time (S2):	12	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.02	kN		
	Max ;	0.07	kN	Max Punch ;	9
Lower Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Ejection Force :	Average ;	0.17	kN		
	Max ;	0.30	kN	Max Punch ;	45

Tablet Data (N/A)

Total Counter :	150086	Tablets			
A.W.C Counter :	132289	Tablets			
HEP :	555 (0.4%)	Tablets	LEP :	232 (0.2%)	Tablets
GOOD :	131502 (99.4%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 20 MG,STKQ1A00
User ID :	rp98927	Batch NO. :	KQM2501920
Print Interval :	50000 Tablets	Running Time :	00 Hour 39 Min 27 Sec
Total Counter :	150086	Total Running Time :	1578 Hour 31 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	15	%	Manual ;	20	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		5500000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		25	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		30.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	12.84	kN	No		
HEP :	20.0 %	11.85	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	10.07	kN	Yes 15	Times	
Ref :		9.88	kN			
LCP :	2.0 %	9.68	kN	Yes 15	Times	
LEP :	20.0 %	7.90	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	6.91	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	1.54	kN	Min :	1.25	kN	Punch NO.;	16
SD :	9.30	%	Max :	1.89	kN	Punch NO.;	55
Mean Main Pressure :	9.84	kN	Min :	8.59	kN	Punch NO.;	16
SD :	6.00	%	Max :	11.03	kN	Punch NO.;	17
Filling Depth :	Increase ;	92	Times	Decrease ;	92	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	30	RPM	Capacity :	234000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	4	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.10	mm	Pre Thickness ;	1.10	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.55	mm	Main Thickness ;	0.55	mm
Filling Depth :	6.67	mm	Current Cam :	Cam B	
Main Air pressure :	0.718	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	21	Min			
Lower Punch Remain Time (S2):	12	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.02	kN		
	Max ;	0.07	kN	Max Punch ;	9
Lower Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Ejection Force :	Average ;	0.17	kN		
	Max ;	0.28	kN	Max Punch ;	45

Tablet Data (N/A)

Total Counter :	150086	Tablets			
A.W.C Counter :	132289	Tablets			
HEP :	555 (0.4%)	Tablets	LEP :	232 (0.2%)	Tablets
GOOD :	131502 (99.4%)	Tablets			