

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 20<br>MG,STKQ1A00 |
| User ID :        | pb97604       | Batch NO. :          | KQM2502120              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 23 Min 21 Sec   |
| Total Counter :  | 5044992       | Total Running Time : | 1719 Hour 25 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 15                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 5500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 40                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 17.05             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 16.39             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 13.37             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 13.11             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 12.85             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 9.83              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 9.18              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |                |       |                |       |    |            |    |
|----------------------|----------------|-------|----------------|-------|----|------------|----|
| Mean Pre Pressure :  | 1.74           | kN    | Min :          | 1.42  | kN | Punch NO.; | 38 |
| SD :                 | 8.40           | %     | Max :          | 2.04  | kN | Punch NO.; | 53 |
| Mean Main Pressure : | 13.28          | kN    | Min :          | 11.29 | kN | Punch NO.; | 38 |
| SD :                 | 6.50           | %     | Max :          | 15.38 | kN | Punch NO.; | 57 |
| Filling Depth :      | Increase ; 269 | Times | Decrease ; 288 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 4        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.36    | mm  | Pre Thickness ;  | 1.37     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.60    | mm  | Main Thickness ; | 0.61     | mm           |
| Filling Depth :               | 6.37    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.679   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 5       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 2       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 14 |
| Lower Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.33 | kN | Max Punch ; | 38 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 55 |

### Tablet Data (N/A)

|                 |                 |         |       |             |         |
|-----------------|-----------------|---------|-------|-------------|---------|
| Total Counter : | 5044992         | Tablets |       |             |         |
| A.W.C Counter : | 5024712         | Tablets |       |             |         |
| HEP :           | 1858 (0.0%)     | Tablets | LEP : | 1540 (0.0%) | Tablets |
| GOOD :          | 5021314 (99.9%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 20<br>MG,STKQ1A00 |
| User ID :        | pb97604       | Batch NO. :          | KQM2502120              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 23 Min 21 Sec   |
| Total Counter :  | 5044992       | Total Running Time : | 1719 Hour 25 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 15                | %       | Manual ;           | 20    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 5500000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 40                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 15                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 30.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 17.77             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 17.09             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 13.94             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 13.67             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 13.40             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 10.25             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 9.57              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.95       | kN  | Min : | 1.67       | kN  | Punch NO.; | 26 |
| SD :                 | 6.10       | %   | Max : | 2.21       | kN  | Punch NO.; | 62 |
| Mean Main Pressure : | 13.74      | kN  | Min : | 11.76      | kN  | Punch NO.; | 26 |
| SD :                 | 5.20       | %   | Max : | 15.20      | kN  | Punch NO.; | 65 |
| Filling Depth :      | Increase ; | 466 | Times | Decrease ; | 411 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 4        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.36    | mm  | Pre Thickness ;  | 1.37     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.66    | mm  | Main Thickness ; | 0.67     | mm           |
| Filling Depth :               | 6.73    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.679   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 5       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 2       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 14 |
| Lower Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.33 | kN | Max Punch ; | 38 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 35 |

### Tablet Data (N/A)

|                 |                 |         |       |             |         |
|-----------------|-----------------|---------|-------|-------------|---------|
| Total Counter : | 5044992         | Tablets |       |             |         |
| A.W.C Counter : | 5024712         | Tablets |       |             |         |
| HEP :           | 1858 (0.0%)     | Tablets | LEP : | 1540 (0.0%) | Tablets |
| GOOD :          | 5021314 (99.9%) | Tablets |       |             |         |