

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 20<br>MG,STKQ1A00 |
| User ID :        | rp98927       | Batch NO. :          | KQM2502220              |
| Print Interval : | 50000 Tablets | Running Time :       | 23 Hour 31 Min 54 Sec   |
| Total Counter :  | 4602424       | Total Running Time : | 1745 Hour 3 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 20                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 5500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 40                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 17.01             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 16.36             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 13.35             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 13.08             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 12.82             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 9.81              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 9.16              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.66  | kN    | Min :      | 1.33  | kN    | Punch NO.; | 26 |
|                      |            | 8.60  | %     | Max :      | 1.96  | kN    | Punch NO.; | 39 |
| Mean Main Pressure : |            | 13.14 | kN    | Min :      | 11.28 | kN    | Punch NO.; | 26 |
|                      | SD :       | 5.60  | %     | Max :      | 15.08 | kN    | Punch NO.; | 65 |
| Filling Depth :      | Increase ; | 1455  | Times | Decrease ; | 1490  | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 5        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.34    | mm  | Pre Thickness ;  | 1.35     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.57    | mm  | Main Thickness ; | 0.58     | mm           |
| Filling Depth :               | 6.05    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.783   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 39      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 9       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.20 | kN | Max Punch ; | 36 |
| Lower Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 38 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 55 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 4602424             | Tablets |       |            |         |
| A.W.C Counter : | 4590849             | Tablets |       |            |         |
| HEP :           | 530 (0.0%)          | Tablets | LEP : | 576 (0.0%) | Tablets |
| GOOD :          | 4589743<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 20<br>MG,STKQ1A00 |
| User ID :        | rp98927       | Batch NO. :          | KQM2502220              |
| Print Interval : | 50000 Tablets | Running Time :       | 23 Hour 31 Min 54 Sec   |
| Total Counter :  | 4602424       | Total Running Time : | 1745 Hour 3 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 20                | %       | Manual ;           | 20           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 5500000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 40                | Min     | Supply Time ;      | 1.0          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 15                | Min     | Supply Time ;      | 2.0          | Sec        |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 17.16             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 16.50             | kN      | Yes                | 10 Rotations | 15 Tablets |
| HCP :                              | 2.0 %       | 13.47             | kN      | Yes                | 15 Times     |            |
| Ref :                              |             | 13.20             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 12.94             | kN      | Yes                | 15 Times     |            |
| LEP :                              | 25.0 %      | 9.90              | kN      | Yes                | 10 Rotations | 15 Tablets |
| LSP :                              | 30.0 %      | 9.24              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |      |       |            |     |            |    |
|----------------------|------------|------|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.81       | kN   | Min : | 1.59       | kN  | Punch NO.; | 20 |
| SD :                 | 6.30       | %    | Max : | 2.10       | kN  | Punch NO.; | 42 |
| Mean Main Pressure : | 13.22      | kN   | Min : | 11.53      | kN  | Punch NO.; | 26 |
| SD :                 | 5.10       | %    | Max : | 15.16      | kN  | Punch NO.; | 42 |
| Filling Depth :      | Increase ; | 1040 | Times | Decrease ; | 988 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 5        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.34    | mm  | Pre Thickness ;  | 1.35     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.63    | mm  | Main Thickness ; | 0.64     | mm           |
| Filling Depth :               | 6.44    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.783   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 39      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 9       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.20 | kN | Max Punch ; | 39 |
| Lower Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 38 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 55 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 4602424             | Tablets |       |            |         |
| A.W.C Counter : | 4590849             | Tablets |       |            |         |
| HEP :           | 530 (0.0%)          | Tablets | LEP : | 576 (0.0%) | Tablets |
| GOOD :          | 4589743<br>(100.0%) | Tablets |       |            |         |