

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Amio 200mg STAI3A00
User ID :	vg96371	Batch NO. :	AIM2500420
Print Interval :	50000 Tablets	Running Time :	01 Hour 48 Min 25 Sec
Total Counter :	139078	Total Running Time :	1751 Hour 5 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	122	%	Manual ;	20	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
		4.5-10mm			Max :	8.00 mm
Target Quantity :		1500000	Tablets			
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		60	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	26.14	kN	No		
HEP :	25.0 %	25.13	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	20.51	kN	Yes	15	Times
Ref :		20.11	kN			
LCP :	2.0 %	19.70	kN	Yes	15	Times
LEP :	25.0 %	15.08	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	14.07	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.17	kN	Min :	2.63	kN	Punch NO.;	50
SD :	9.20	%	Max :	3.92	kN	Punch NO.;	28
Mean Main Pressure :	19.91	kN	Min :	17.81	kN	Punch NO.;	10
SD :	4.70	%	Max :	21.97	kN	Punch NO.;	22
Filling Depth :	Increase ;	103	Times	Decrease ;	112	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	18	RPM
Pre Upper Punch Position ;	2.20	mm			
Pre Lower Punch Position ;	3.40	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	2.80	mm			
Main Lower Punch Position ;	3.40	mm	Main Thickness ;	0.60	mm
Filling Depth :	6.69	mm	Current Cam :	Cam B	
Main Air pressure :	0.807	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	36	Min			
Lower Punch Remain Time (S2):	7	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.30	kN		
	Max ;	0.38	kN	Max Punch ;	49
Lower Punch :	Average ;	0.10	kN		
	Max ;	0.16	kN	Max Punch ;	35
Ejection Force :	Average ;	0.03	kN		
	Max ;	0.04	kN	Max Punch ;	2

Tablet Data (N/A)

Total Counter :	139078	Tablets			
A.W.C Counter :	143144	Tablets			
HEP :	263 (0.2%)	Tablets	LEP :	126 (0.1%)	Tablets
GOOD :	142755 (99.7%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Amio 200mg STAI3A00
User ID :	vg96371	Batch NO. :	AIM2500420
Print Interval :	50000 Tablets	Running Time :	01 Hour 48 Min 25 Sec
Total Counter :	139078	Total Running Time :	1751 Hour 5 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	124	%	Manual ;	20	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
		4.5-10mm			Max :	8.00 mm
Target Quantity :		1500000	Tablets			
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		60	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	26.47	kN	No		
HEP :	25.0 %	25.45	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	20.77	kN	Yes 15	Times	
Ref :		20.36	kN			
LCP :	2.0 %	19.96	kN	Yes 15	Times	
LEP :	25.0 %	15.27	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	14.25	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.58	kN	Min :	1.87	kN	Punch NO.;	13
SD :	11.00	%	Max :	3.27	kN	Punch NO.;	45
Mean Main Pressure :	19.86	kN	Min :	16.99	kN	Punch NO.;	13
SD :	5.10	%	Max :	22.68	kN	Punch NO.;	45
Filling Depth :	Increase ; 94	Times	Decrease ;	100	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	19	RPM
Pre Upper Punch Position ;	2.20	mm			
Pre Lower Punch Position ;	3.40	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	2.80	mm			
Main Lower Punch Position ;	3.40	mm	Main Thickness ;	0.60	mm
Filling Depth :	6.85	mm	Current Cam :	Cam B	
Main Air pressure :	0.807	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	36	Min			
Lower Punch Remain Time (S2):	7	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.30	kN		
	Max ;	0.43	kN	Max Punch ;	49
Lower Punch :	Average ;	0.10	kN		
	Max ;	0.16	kN	Max Punch ;	35
Ejection Force :	Average ;	0.03	kN		
	Max ;	0.04	kN	Max Punch ;	2

Tablet Data (N/A)

Total Counter :	139078	Tablets			
A.W.C Counter :	143144	Tablets			
HEP :	263 (0.2%)	Tablets	LEP :	126 (0.1%)	Tablets
GOOD :	142755 (99.7%)	Tablets			