

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Amio 200mg STAI3A00   |
| User ID :        | pb97604       | Batch NO. :          | AIM2500420            |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 52 Min 37 Sec |
| Total Counter :  | 909882        | Total Running Time : | 1759 Hour 9 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |       |                      |
|------------------------------------|-------------|----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ;      | 120      | %       | Manual ;           | 20    | RPM                  |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min : | 2.50 mm              |
|                                    |             | 4.5-10mm |         |                    | Max : | 8.00 mm              |
| Target Quantity :                  |             | 1500000  | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |       |                      |
| Oil Lubrication :                  |             |          |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |             | 60       | Min     | Supply Time ;      | 1.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |       |                      |
| Initial Reject Time :              |             | 2        | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |       |                      |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |       |                      |
|                                    |             |          |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 %      | 27.62    | kN      | No                 |       |                      |
| HEP :                              | 25.0 %      | 26.56    | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %       | 21.67    | kN      | Yes                | 15    | Times                |
| Ref :                              |             | 21.25    | kN      |                    |       |                      |
| LCP :                              | 2.0 %       | 20.82    | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 %      | 15.94    | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 %      | 14.87    | kN      | No                 |       |                      |
| SD limit :                         | 35.0 %      |          |         | Yes                |       |                      |
| Pre HSP                            |             | 30.00    | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None        |          |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.46       | kN  | Min : | 2.73       | kN  | Punch NO.; | 14 |
| SD :                 | 11.70      | %   | Max : | 4.63       | kN  | Punch NO.; | 2  |
| Mean Main Pressure : | 21.34      | kN  | Min : | 18.62      | kN  | Punch NO.; | 14 |
| SD :                 | 6.30       | %   | Max : | 25.25      | kN  | Punch NO.; | 2  |
| Filling Depth :      | Increase ; | 581 | Times | Decrease ; | 595 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 2.20    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 3.40    | mm  | Pre Thickness ;  | 1.20     | mm           |
| Main Upper Punch Position ;   | 2.80    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.42    | mm  | Main Thickness ; | 0.62     | mm           |
| Filling Depth :               | 6.42    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.816   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 35      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 6       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.09 | kN |             |   |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 1 |
| Lower Punch :    | Average ; | 0.26 | kN |             |   |
|                  | Max ;     | 0.74 | kN | Max Punch ; | 1 |
| Ejection Force : | Average ; | 0.03 | kN |             |   |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 3 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 909882            | Tablets |       |            |         |
| A.W.C Counter : | 906996            | Tablets |       |            |         |
| HEP :           | 558 (0.1%)        | Tablets | LEP : | 218 (0.0%) | Tablets |
| GOOD :          | 906220<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Amio 200mg STAI3A00   |
| User ID :        | pb97604       | Batch NO. :          | AIM2500420            |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 52 Min 37 Sec |
| Total Counter :  | 909882        | Total Running Time : | 1759 Hour 9 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 124      | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 29.32    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 28.20    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 23.01    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 22.56    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 22.11    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 16.92    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 15.79    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.16       | kN  | Min : | 2.39       | kN  | Punch NO.; | 33 |
| SD :                 | 13.60      | %   | Max : | 4.34       | kN  | Punch NO.; | 2  |
| Mean Main Pressure : | 22.17      | kN  | Min : | 19.78      | kN  | Punch NO.; | 12 |
| SD :                 | 5.60       | %   | Max : | 25.38      | kN  | Punch NO.; | 2  |
| Filling Depth :      | Increase ; | 271 | Times | Decrease ; | 290 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 19       | RPM          |
| Pre Upper Punch Position ;    | 2.20    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 3.40    | mm  | Pre Thickness ;  | 1.20     | mm           |
| Main Upper Punch Position ;   | 2.80    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.42    | mm  | Main Thickness ; | 0.62     | mm           |
| Filling Depth :               | 6.79    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.816   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 35      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 6       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.09 | kN |             |   |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 1 |
| Lower Punch :    | Average ; | 0.26 | kN |             |   |
|                  | Max ;     | 0.55 | kN | Max Punch ; | 1 |
| Ejection Force : | Average ; | 0.03 | kN |             |   |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 3 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 909882            | Tablets |       |            |         |
| A.W.C Counter : | 906996            | Tablets |       |            |         |
| HEP :           | 558 (0.1%)        | Tablets | LEP : | 218 (0.0%) | Tablets |
| GOOD :          | 906220<br>(99.9%) | Tablets |       |            |         |