

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | GLICLAZIDE MR 30MG    |
| User ID :        | rp98927       | Batch NO. :          | GLC TAB XV 25 003     |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 01 Min 13 Sec |
| Total Counter :  | 757158        | Total Running Time : | 1769 Hour 40 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 78       | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 14.12    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 13.58    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 11.08    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 10.86    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 10.64    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 8.15     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 7.60     | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.37  | kN    | Min :      | 1.27  | kN    | Punch NO.; | 9  |
|                      |            | 3.50  | %     | Max :      | 1.50  | kN    | Punch NO.; | 46 |
| Mean Main Pressure : |            | 10.93 | kN    | Min :      | 10.02 | kN    | Punch NO.; | 9  |
|                      | SD :       | 3.10  | %     | Max :      | 11.67 | kN    | Punch NO.; | 46 |
| Filling Depth :      | Increase ; | 23    | Times | Decrease ; | 59    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 16       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.32    | mm  | Pre Thickness ;  | 2.33     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.48    | mm  | Main Thickness ; | 1.49     | mm           |
| Filling Depth :               | 5.33    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.763   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 52      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 23      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.02 | kN | Max Punch ; | 28 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 13 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 7  |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 757158             | Tablets |       |        |         |
| A.W.C Counter : | 746162             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 746162<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | GLICLAZIDE MR 30MG    |
| User ID :        | rp98927       | Batch NO. :          | GLC TAB XV 25 003     |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 01 Min 13 Sec |
| Total Counter :  | 757158        | Total Running Time : | 1769 Hour 40 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 78       | %       | Manual ;           | 24           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 13.74    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 13.21    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 10.78    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 10.57    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 10.36    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 7.93     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 7.40     | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.18  | kN    | Min :      | 1.01  | kN    | Punch NO.; | 18 |
|                      |            | 6.40  | %     | Max :      | 1.38  | kN    | Punch NO.; | 7  |
| Mean Main Pressure : |            | 10.42 | kN    | Min :      | 9.49  | kN    | Punch NO.; | 18 |
|                      | SD :       | 4.20  | %     | Max :      | 11.38 | kN    | Punch NO.; | 6  |
| Filling Depth :      | Increase ; | 39    | Times | Decrease ; | 48    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 16       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.40    | mm  | Pre Thickness ;  | 2.41     | mm           |
| Main Upper Punch Position ;   | 2.98    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.56    | mm  | Main Thickness ; | 1.58     | mm           |
| Filling Depth :               | 5.56    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.763   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 52      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 23      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 36 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 13 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 7  |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 757158             | Tablets |       |        |         |
| A.W.C Counter : | 746162             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 746162<br>(100.0%) | Tablets |       |        |         |