

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40 MG<br>STKQ1B00 |
| User ID :        | os190259      | Batch NO. :          | KQM2507440              |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 11 Min 29 Sec   |
| Total Counter :  | 711978        | Total Running Time : | 1802 Hour 21 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 55                | %       | Manual ;           | 10    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 2750000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 35                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 35.50             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 34.14             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 27.86             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 27.31             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 26.76             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 20.48             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 19.12             | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.17  | kN    | Min :      | 1.11  | kN    | Punch NO.; | 38 |
|                      |            | 2.40  | %     | Max :      | 1.23  | kN    | Punch NO.; | 43 |
| Mean Main Pressure : |            | 27.40 | kN    | Min :      | 24.94 | kN    | Punch NO.; | 28 |
|                      | SD :       | 4.40  | %     | Max :      | 30.20 | kN    | Punch NO.; | 43 |
| Filling Depth :      | Increase ; | 51    | Times | Decrease ; | 48    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.00    | mm  | Pre Thickness ;  | 2.01     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.94    | mm  | Main Thickness ; | 0.94     | mm           |
| Filling Depth :               | 4.68    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.801   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 24      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 10      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.28 | kN | Max Punch ; | 24 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 26 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 34 |

### Tablet Data (N/A)

|                 |                |         |       |            |         |
|-----------------|----------------|---------|-------|------------|---------|
| Total Counter : | 711978         | Tablets |       |            |         |
| A.W.C Counter : | 685990         | Tablets |       |            |         |
| HEP :           | 157 (0.0%)     | Tablets | LEP : | 521 (0.1%) | Tablets |
| GOOD :          | 685312 (99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40 MG<br>STKQ1B00 |
| User ID :        | os190259      | Batch NO. :          | KQM2507440              |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 11 Min 29 Sec   |
| Total Counter :  | 711978        | Total Running Time : | 1802 Hour 21 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 55                | %       | Manual ;           | 10           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 35                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 35.67             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 34.29             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 27.98             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 27.44             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 26.89             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 20.58             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 19.20             | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.13       | kN | Min : | 1.09       | kN | Punch NO.; | 12 |
| SD :                 | 1.80       | %  | Max : | 1.18       | kN | Punch NO.; | 9  |
| Mean Main Pressure : | 27.56      | kN | Min : | 25.59      | kN | Punch NO.; | 63 |
| SD :                 | 2.90       | %  | Max : | 29.54      | kN | Punch NO.; | 26 |
| Filling Depth :      | Increase ; | 14 | Times | Decrease ; | 32 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.00    | mm  | Pre Thickness ;  | 2.00     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.10    | mm  | Main Thickness ; | 1.11     | mm           |
| Filling Depth :               | 4.32    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.801   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 24      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 10      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.30 | kN | Max Punch ; | 24 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 26 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 34 |

### Tablet Data (N/A)

|                 |                |         |       |            |         |
|-----------------|----------------|---------|-------|------------|---------|
| Total Counter : | 711978         | Tablets |       |            |         |
| A.W.C Counter : | 685990         | Tablets |       |            |         |
| HEP :           | 157 (0.0%)     | Tablets | LEP : | 521 (0.1%) | Tablets |
| GOOD :          | 685312 (99.9%) | Tablets |       |            |         |