

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA 50MG<br>STSH4A00 |
| User ID :        | pb97604       | Batch NO. :          | SHM2502450              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 52 Min 06 Sec   |
| Total Counter :  | 1000824       | Total Running Time : | 1961 Hour 31 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 46                | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 6000000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 16.78             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 16.14             | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 1.5 %       | 13.10             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 12.91             | kN      |                    |              |           |
| LCP :                              | 1.5 %       | 12.71             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 9.68              | kN      | No 35              | Rotations 35 | Tablets   |
| LSP :                              | 30.0 %      | 9.04              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.97       | kN  | Min : | 1.88       | kN  | Punch NO.; | 9  |
| SD :                 | 1.70       | %   | Max : | 2.04       | kN  | Punch NO.; | 30 |
| Mean Main Pressure : | 12.97      | kN  | Min : | 12.06      | kN  | Punch NO.; | 19 |
| SD :                 | 3.10       | %   | Max : | 13.84      | kN  | Punch NO.; | 57 |
| Filling Depth :      | Increase ; | 304 | Times | Decrease ; | 384 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 13                | RPM          |
| Pre Upper Punch Position ;    | 3.00   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.60   | mm  | Pre Thickness ;  | 2.60              | mm           |
| Main Upper Punch Position ;   | 3.00   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.13   | mm  | Main Thickness ; | 2.13              | mm           |
| Filling Depth :               | 5.95   | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.786  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 57     | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 16     | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 63 |
| Lower Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 70 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1000824            | Tablets |       |            |         |
| A.W.C Counter : | 995745             | Tablets |       |            |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 198 (0.0%) | Tablets |
| GOOD :          | 995547<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA 50MG<br>STSH4A00 |
| User ID :        | pb97604       | Batch NO. :          | SHM2502450              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 52 Min 06 Sec   |
| Total Counter :  | 1000824       | Total Running Time : | 1961 Hour 31 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 42                | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 6000000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 16.23             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 15.60             | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 1.5 %  | 12.67             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 12.48             | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 12.30             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 9.36              | kN      | Yes 35             | Rotations 35 | Tablets   |
| LSP :                              | 30.0 % | 8.74              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.96       | kN  | Min : | 1.82       | kN  | Punch NO.; | 20 |
| SD :                 | 2.80       | %   | Max : | 2.09       | kN  | Punch NO.; | 74 |
| Mean Main Pressure : | 12.45      | kN  | Min : | 11.17      | kN  | Punch NO.; | 20 |
| SD :                 | 4.20       | %   | Max : | 13.63      | kN  | Punch NO.; | 50 |
| Filling Depth :      | Increase ; | 310 | Times | Decrease ; | 387 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 3.00   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.50              | mm           |
| Main Upper Punch Position ;   | 3.00   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.23   | mm  | Main Thickness ; | 2.23              | mm           |
| Filling Depth :               | 5.60   | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.786  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 57     | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 16     | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 63 |
| Lower Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 70 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1000824            | Tablets |       |            |         |
| A.W.C Counter : | 995745             | Tablets |       |            |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 198 (0.0%) | Tablets |
| GOOD :          | 995547<br>(100.0%) | Tablets |       |            |         |