

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	METHO 500MG STHA1A00
User ID :	rp98927	Batch NO. :	HAM2501050
Print Interval :	50000 Tablets	Running Time :	02 Hour 37 Min 58 Sec
Total Counter :	250070	Total Running Time :	2076 Hour 37 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	98	%	Manual ;	30	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		750000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		80	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	13.39	kN	No		
HEP :	25.0 %	12.88	kN	Yes 10	Rotations 20	Tablets
HCP :	1.5 %	10.46	kN	Yes 15	Times	
Ref :		10.30	kN			
LCP :	1.5 %	10.15	kN	Yes 15	Times	
LEP :	25.0 %	7.73	kN	No 35	Rotations 35	Tablets
LSP :	30.0 %	7.21	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	1 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	0.99	kN	Min :	0.96	kN	Punch NO.;	20
		0.90	%	Max :	1.01	kN	Punch NO.;	22
Mean Main Pressure :		10.15	kN	Min :	8.21	kN	Punch NO.;	10
	SD :	7.20	%	Max :	11.69	kN	Punch NO.;	43
Filling Depth :	Increase ;	436	Times	Decrease ;	377	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	15	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	6.50	mm	Pre Thickness ;	4.50	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	4.50	mm	Main Thickness ;	2.50	mm
Filling Depth :	8.93	mm	Current Cam :	Cam C	
Main Air pressure :	0.800	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	2	Min			
Lower Punch Remain Time (S2):	24	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.06	kN		
	Max ;	0.09	kN	Max Punch ;	40
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.13	kN	Max Punch ;	47
Ejection Force :	Average ;	0.11	kN		
	Max ;	0.17	kN	Max Punch ;	44

Tablet Data (N/A)

Total Counter :	250070	Tablets			
A.W.C Counter :	224858	Tablets			
HEP :	377 (0.2%)	Tablets	LEP :	670 (0.3%)	Tablets
GOOD :	223811 (99.5%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	METHO 500MG STHA1A00
User ID :	rp98927	Batch NO. :	HAM2501050
Print Interval :	50000 Tablets	Running Time :	02 Hour 37 Min 58 Sec
Total Counter :	250070	Total Running Time :	2076 Hour 37 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	94	%	Manual ;	30	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		750000	Tablets			
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		80	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	12.44	kN	No		
HEP :	25.0 %	11.96	kN	Yes 10	Rotations 20	Tablets
HCP :	1.5 %	9.71	kN	Yes 15	Times	
Ref :		9.57	kN			
LCP :	1.5 %	9.43	kN	Yes 15	Times	
LEP :	25.0 %	7.18	kN	Yes 35	Rotations 35	Tablets
LSP :	30.0 %	6.70	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	0.75	kN	Min :	0.73	kN	Punch NO.;	2
SD :	2.20	%	Max :	0.80	kN	Punch NO.;	13
Mean Main Pressure :	9.29	kN	Min :	8.31	kN	Punch NO.;	7
SD :	4.70	%	Max :	10.50	kN	Punch NO.;	35
Filling Depth :	Increase ;	356	Times	Decrease ;	352	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	6.50	mm	Pre Thickness ;	4.50	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	4.66	mm	Main Thickness ;	2.66	mm
Filling Depth :	8.14	mm	Current Cam :	Cam C	
Main Air pressure :	0.800	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	2	Min			
Lower Punch Remain Time (S2):	24	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.06	kN		
	Max ;	0.09	kN	Max Punch ;	40
Lower Punch :	Average ;	0.09	kN		
	Max ;	0.13	kN	Max Punch ;	47
Ejection Force :	Average ;	0.11	kN		
	Max ;	0.18	kN	Max Punch ;	44

Tablet Data (N/A)

Total Counter :	250070	Tablets			
A.W.C Counter :	224858	Tablets			
HEP :	377 (0.2%)	Tablets	LEP :	670 (0.3%)	Tablets
GOOD :	223811 (99.5%)	Tablets			