

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU 800MG STIU2D00
User ID :	sj97904	Batch NO. :	IUM2500980
Print Interval :	50000 Tablets	Running Time :	00 Hour 31 Min 52 Sec
Total Counter :	50008	Total Running Time :	2124 Hour 0 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	170	%	Manual ;	30	RPM
Filling Cam :		Cam D 12.5- 18mm		Cam Range :	Min :	10.50 mm
Target Quantity :		150000	Tablets		Max :	16.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	28.24	kN	No		
HEP :	25.0 %	27.16	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	22.16	kN	Yes 15	Times	
Ref :		21.73	kN			
LCP :	2.0 %	21.29	kN	Yes 15	Times	
LEP :	25.0 %	16.29	kN	No 15	Rotations 20	Tablets
LSP :	30.0 %	15.21	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.41	kN	Min :	3.29	kN	Punch NO.;	14
SD :	1.50	%	Max :	3.54	kN	Punch NO.;	42
Mean Main Pressure :	21.16	kN	Min :	19.95	kN	Punch NO.;	41
SD :	2.60	%	Max :	22.07	kN	Punch NO.;	18
Filling Depth :	Increase ; 75	Times	Decrease ; 57	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	26	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	9.00	mm	Pre Thickness ;	6.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	7.70	mm	Main Thickness ;	4.71	mm
Filling Depth :	11.89	mm	Current Cam :	Cam D 12.5- 18mm	
Main Air pressure :	0.712	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	9	Min			
Lower Punch Remain Time (S2):	5	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.42	kN		
	Max ;	0.69	kN	Max Punch ;	5
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.17	kN	Max Punch ;	19
Ejection Force :	Average ;	1.19	kN		
	Max ;	1.45	kN	Max Punch ;	17

Tablet Data (N/A)

Total Counter :	50008	Tablets			
A.W.C Counter :	48560	Tablets			
HEP :	2 (0.0%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	48557 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU 800MG STIU2D00
User ID :	sj97904	Batch NO. :	IUM2500980
Print Interval :	50000 Tablets	Running Time :	00 Hour 31 Min 52 Sec
Total Counter :	50008	Total Running Time :	2124 Hour 0 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	174	%	Manual ;	30	RPM
Filling Cam :		Cam D 12.5- 18mm		Cam Range :	Min :	10.50 mm
Target Quantity :		150000	Tablets		Max :	16.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	29.22	kN	No		
HEP :	25.0 %	28.09	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	22.92	kN	Yes 20	Times	
Ref :		22.47	kN			
LCP :	2.0 %	22.02	kN	Yes 20	Times	
LEP :	25.0 %	16.86	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	15.73	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.78	kN	Min :	3.67	kN	Punch NO.;	9
SD :	1.60	%	Max :	3.92	kN	Punch NO.;	48
Mean Main Pressure :	24.35	kN	Min :	22.75	kN	Punch NO.;	30
SD :	3.40	%	Max :	27.05	kN	Punch NO.;	51
Filling Depth :	Increase ;	147	Times	Decrease ;	145	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	26	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	9.00	mm	Pre Thickness ;	6.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	7.85	mm	Main Thickness ;	4.86	mm
Filling Depth :	11.88	mm	Current Cam :	Cam D 12.5- 18mm	
Main Air pressure :	0.712	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	9	Min			
Lower Punch Remain Time (S2):	5	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.40	kN		
	Max ;	0.63	kN	Max Punch ;	5
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.17	kN	Max Punch ;	19
Ejection Force :	Average ;	1.19	kN		
	Max ;	1.45	kN	Max Punch ;	17

Tablet Data (N/A)

Total Counter :	50008	Tablets			
A.W.C Counter :	48560	Tablets			
HEP :	2 (0.0%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	48557 (100.0%)	Tablets			