

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MET 1000mg STYQ1A00
User ID :	pb97604	Batch NO. :	YQM2502110
Print Interval :	50000 Tablets	Running Time :	01 Hour 03 Min 38 Sec
Total Counter :	100016	Total Running Time :	2147 Hour 53 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	135	%	Manual ;	30	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2250000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	27.96	kN	No		
HEP :	25.0 %	26.88	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	21.94	kN	Yes 15	Times	
Ref :		21.51	kN			
LCP :	2.0 %	21.08	kN	Yes 15	Times	
LEP :	25.0 %	16.13	kN	No 15	Rotations 20	Tablets
LSP :	30.0 %	15.05	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	3.35	kN	Min :	3.13	kN	Punch NO.;	17
		2.00	%	Max :	3.48	kN	Punch NO.;	12
Mean Main Pressure :		21.87	kN	Min :	20.67	kN	Punch NO.;	18
	SD :	2.40	%	Max :	23.07	kN	Punch NO.;	12
Filling Depth :	Increase ;	20	Times	Decrease ;	32	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	20	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	8.30	mm	Pre Thickness ;	5.31	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	7.33	mm	Main Thickness ;	4.34	mm
Filling Depth :	10.48	mm	Current Cam :	Cam C	
Main Air pressure :	0.731	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	8	Min			
Lower Punch Remain Time (S2):	5	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.08	kN		
	Max ;	0.18	kN	Max Punch ;	25
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.11	kN	Max Punch ;	32
Ejection Force :	Average ;	0.03	kN		
	Max ;	0.11	kN	Max Punch ;	1

Tablet Data (N/A)

Total Counter :	100016	Tablets			
A.W.C Counter :	97763	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	97763 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MET 1000mg STYQ1A00
User ID :	pb97604	Batch NO. :	YQM2502110
Print Interval :	50000 Tablets	Running Time :	01 Hour 03 Min 38 Sec
Total Counter :	100016	Total Running Time :	2147 Hour 53 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	135	%	Manual ;	30	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2250000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	27.59	kN	No		
HEP :	25.0 %	26.52	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	21.64	kN	Yes 20	Times	
Ref :		21.22	kN			
LCP :	2.0 %	20.79	kN	Yes 20	Times	
LEP :	25.0 %	15.91	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	14.85	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.47	kN	Min :	3.30	kN	Punch NO.;	17
SD :	2.10	%	Max :	3.69	kN	Punch NO.;	41
Mean Main Pressure :	21.81	kN	Min :	20.23	kN	Punch NO.;	17
SD :	3.60	%	Max :	24.09	kN	Punch NO.;	41
Filling Depth :	Increase ; 29	Times	Decrease ; 21	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	20	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	8.30	mm	Pre Thickness ;	5.31	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	7.46	mm	Main Thickness ;	4.47	mm
Filling Depth :	10.14	mm	Current Cam :	Cam C	
Main Air pressure :	0.731	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	8	Min			
Lower Punch Remain Time (S2):	5	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.08	kN		
	Max ;	0.16	kN	Max Punch ;	27
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.11	kN	Max Punch ;	32
Ejection Force :	Average ;	0.03	kN		
	Max ;	0.11	kN	Max Punch ;	1

Tablet Data (N/A)

Total Counter :	100016	Tablets			
A.W.C Counter :	97763	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	97763 (100.0%)	Tablets			