

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MIRTAZAPINE 15 MG     |
| User ID :        | ks99394       | Batch NO. :          | MRT TAB XV 25 001B    |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 48 Min 00 Sec |
| Total Counter :  | 150086        | Total Running Time : | 2208 Hour 2 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 20       | %       | Manual ;           | 57           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 513333   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 70       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 12.20    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 11.73    | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 9.57     | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 9.38     | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 9.19     | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 7.04     | kN      | No 15              | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 6.57     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |              |       |            |      |       |            |    |
|----------------------|--------------|-------|------------|------|-------|------------|----|
| Mean Pre Pressure :  | 2.42         | kN    | Min :      | 2.36 | kN    | Punch NO.; | 40 |
| SD :                 | 1.00         | %     | Max :      | 2.49 | kN    | Punch NO.; | 77 |
| Mean Main Pressure : | 9.34         | kN    | Min :      | 9.15 | kN    | Punch NO.; | 56 |
| SD :                 | 0.90         | %     | Max :      | 9.55 | kN    | Punch NO.; | 48 |
| Filling Depth :      | Increase ; 4 | Times | Decrease ; | 12   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 231000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 5        | RPM          |
| Pre Upper Punch Position ;    | 4.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.65    | mm  | Pre Thickness ;  | 2.15     | mm           |
| Main Upper Punch Position ;   | 4.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.35    | mm  | Main Thickness ; | 1.85     | mm           |
| Filling Depth :               | 7.39    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.750   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 49      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.22 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 59 |
| Ejection Force : | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 56 |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 150086            | Tablets |       |        |         |
| A.W.C Counter : | 129684            | Tablets |       |        |         |
| HEP :           | 376 (0.3%)        | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 129308<br>(99.7%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MIRTAZAPINE 15 MG     |
| User ID :        | ks99394       | Batch NO. :          | MRT TAB XV 25 001B    |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 48 Min 00 Sec |
| Total Counter :  | 150086        | Total Running Time : | 2208 Hour 2 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 20       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 513333   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 70       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 12.09    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 11.63    | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 9.49     | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 9.30     | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 9.12     | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 %      | 6.98     | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 6.51     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |              |       |            |      |       |            |    |
|----------------------|--------------|-------|------------|------|-------|------------|----|
| Mean Pre Pressure :  | 2.43         | kN    | Min :      | 2.37 | kN    | Punch NO.; | 23 |
| SD :                 | 1.30         | %     | Max :      | 2.51 | kN    | Punch NO.; | 11 |
| Mean Main Pressure : | 9.09         | kN    | Min :      | 8.54 | kN    | Punch NO.; | 21 |
| SD :                 | 2.60         | %     | Max :      | 9.59 | kN    | Punch NO.; | 39 |
| Filling Depth :      | Increase ; 5 | Times | Decrease ; | 20   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 231000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 5        | RPM          |
| Pre Upper Punch Position ;    | 4.51    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.80    | mm  | Pre Thickness ;  | 2.29     | mm           |
| Main Upper Punch Position ;   | 4.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.52    | mm  | Main Thickness ; | 2.02     | mm           |
| Filling Depth :               | 6.40    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.750   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 49      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.22 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 59 |
| Ejection Force : | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 39 |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 150086            | Tablets |       |        |         |
| A.W.C Counter : | 129684            | Tablets |       |        |         |
| HEP :           | 376 (0.3%)        | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 129308<br>(99.7%) | Tablets |       |        |         |